



ITEM # **TNR TN1550**

## FEATURES:

- Left and Right Thermal Chambers, Forceps Wells, Dispenser, Hot Plate, and Cryo Plate are all Individually Temperature Controlled
- Histo Trimmers on Both Sides for Left/Right Handed Users
- Large Adjustable Magnifying Glass
- Automatic and Manual Modes of Operation
- Bright, Push Button Adjustable LED Lighting
- 6 Liter Capacity Paraffin Tank
- Adjustable Temperature Range for Heated Areas
- Large Cryo Module with Adjustable Temperature Range
- Foot Pedal or Hand Trigger to Release Paraffin
- Fully Programmable ON/OFF is User Friendly
- Oversized, 4-Line Display Screen Enables Viewing at Any Angle and Viewing of All Temperature Settings
- Large, Easy-to-Open Drainage Trays Under All Work Areas
- Heated Forceps Holder
- Hinged, Backward Sliding Storage Covers

# TN1550 EMBEDDING STATION

## TECHNICAL SPECIFICATIONS:

- Electric Security Type: Class I
- (Ambient) Temp range: 85°C
- Cryo Plate Temp range : 30°C ~ -15°C
- Temperature tolerance: ±10°C
- Paraffin Tank: 495m x 132 x 95mm
- L/R Thermal Chambers: 235 x 160 x 50mm
- Working Plate: 513 x 93mm
- Voltage: AC110V 60HZ
- Power: 1000VA
- Overall Dimension: 680 x 580 x 420mm
- Weight: 35kg

**Mercedes**  
MEDICAL

PHONE: 800.331.2716 | FAX: 800.359.8807 | WEB: MERCEDESMEDICAL.COM

All offers are subject to availability and are subject to change without notice. The names of various manufacturers, their instruments and their products referred to herein may be protected by trademark or other law, and are used herein solely for purpose of reference. All trademarks are property of their respective owners. Mercedes Medical expressly disclaims any affiliation with or sponsorship referred to within. To stop receiving faxes, please fill out the information below and fax back to the number above. Name \_\_\_\_\_ Fax Number \_\_\_\_\_

58128