

item#
MES CHUCK30

item#
MES CHUCK40

Qty: Each

SPECIMEN DISCS

Sizes: 30mm or 40mm



- Made of aluminum with comfort grips
- Embossed circular discs with “o-ring” for easy handling and edge protection
- For use on Avantik® QS11/QS12 Cryostat and Microm HM520/HM525 Cryostats
- Used in conjunction with heat extractor
- Helps to flatten and freeze samples

ALSO AVAILABLE: LEICA® CRYOSTAT SPECIMEN DISCS

REQUEST SAMPLE TODAY!

Mercedes
MEDICAL

PHONE: 800.331.2716 | FAX: 800.359.8807 | WEB: MERCEDESMEDICAL.COM

All offers are subject to availability and are subject to change without notice. The names of various manufacturers, their instruments and their products referred to herein may be protected by trademark or other law, and are used herein solely for purpose of reference. Mercedes Medical expressly disclaims any affiliation with or sponsorship by the companies referred to herein. All trademarks are property of their respective owners. To stop receiving faxes, please fill out the information below and fax back to the number above. **Name** _____ **Fax Number** _____

62155