

TANNER  
scientific

# path pouch™

**BLUE BIOPSY BAGS**

**THE ONLY  
BLUE BAG  
ON THE  
MARKET!**



**ITEM #  
TNR PATHBAGBLU**

**Size: 30x45mm (Small)**

- Manufactured from solvent resistant nylon
- Thin mesh (0.2 x 0.2 mm) provides excellent fluid exchange and avoids small sample loss
- Suitable for processing small biopsy and histological specimens
- Convenient peel-open bag allows for placing and removing specimen rapidly
- **Packaged 100/Bag, 10 Bags/Box (1,000/Box)**

**Mercedes**  
MEDICAL

PHONE: 800.331.2716 | FAX: 800.359.8807 | WEB: MERCEDESMEDICAL.COM

All offers are subject to availability and are subject to change without notice. The names of various manufacturers, their instruments and their products referred to herein may be protected by trademark or other law, and are used herein solely for purpose of reference. Mercedes Medical expressly disclaims any affiliation with or sponsorship by the companies referred to herein. Tanner Scientific and all associated marks are trademarks of Tanner Scientific LLC. To stop receiving faxes, please fill out the information and fax back to the number above. Name \_\_\_\_\_ Fax Number \_\_\_\_\_

67677  
REV012618