

BLUWING™



Safety Butterflies

w/Adapter & 12" Tubing

item #: MER 1921 - 21G

item #: MER 1923 - 23G

item #: MER 1925 - 25G

50/Box



NEW & IMPROVED

- High quality processing and control
- Reliable safety house for holding the cannula sharp point
- Small safety house design for more convenient use
- Semi-hard safety house is more comfortable for patients
- **Softer, more elastic tubing and sharper needle for easier handling**
- Flexible and kink-resistant elastic tube
- Just one step to lock the needle by drawing back the tube after use, this action can be taken by **one hand**
- Sterilized with ethylene oxide



PHONE: 800.331.2716 | FAX: 800.359.8807 | WEB: MERCEDESMEDICAL.COM