

# BLUWING™



## Safety Butterflies

w/Adapter & 12" Tubing

item #: MER 1921 - 21G

item #: MER 1923 - 23G

item #: MER 1925 - 25G

50/Box



### NEW & IMPROVED

- High quality processing and control
- Reliable safety house for holding the cannula sharp point
- Small safety house design for more convenient use
- Semi-hard safety house is more comfortable for patients
- **Softer, more elastic tubing and sharper needle for easier handling**
- Flexible and kink-resistant elastic tube
- Just one step to lock the needle by drawing back the tube after use, this action can be taken by **one hand**
- Sterilized with ethylene oxide



PHONE: 800.331.2716 | FAX: 800.359.8807 | WEB: MERCEDESMEDICAL.COM



## BD® VACUTAINER™ PLUS PLASTIC BLOOD COLLECTION TUBES

BD Vacutainer® Plus Plastic Tubes offer a safe method for blood collection. The efficacy, performance, and ease of use of BD Vacutainer® Plus Plastic Tubes has been clinically proven and is supported by more than 100 clinical studies, involving a vast array of analytes and diagnostic instruments. Disposal of these plastic tubes is safe, simple and in accordance with EPA guidelines. Plastic tubes reduce the risk of tube breakage and specimen spillage.

ITEM #	TUBE SIZE (MM)	CAP COLOR	DRAW VOLUME	DESCRIPTION	QTY
<b>BEC 366668</b>	13 x 75		3.0 mL	Clot Activator Silicone Coated, Serum	100/Rack
<b>BEC 367812</b>	13 x 75		4.0 mL	Clot Activator Silicone Coated, Serum	100/Rack
<b>BEC 367815</b>	13 x 100		6.0 mL	Clot Activator Silicone Coated, Serum	100/Rack
<b>BEC 367820</b>	16 x 100		10.0 mL	Clot Activator Silicone Coated, Serum	100/Rack
<b>BEC 367842</b>	13 x 75		2.0 mL	Spray-Coated K2EDTA, Whole Blood	100/Rack
<b>BEC 367899</b>	13 x 100		6.0 mL	Spray-Coated K2EDTA, Whole Blood	100/Rack
<b>BEC 367841</b>	13 x 75		2.0 mL	Spray-Coated K2EDTA, Whole Blood	100/Rack
<b>BEC 367844</b>	13 x 75		4.0 mL	Spray-Coated K2EDTA, Whole Blood	100/Rack
<b>BEC 367856</b>	13 x 75		3.0 mL	Spray-Coated K2EDTA, Whole Blood	100/Rack
<b>BEC 367861</b>	13 x 75		4.0 mL	Spray-Coated K2EDTA, Whole Blood	100/Rack
<b>BEC 367863</b>	13 x 100		6.0 mL	Spray-Coated K2EDTA, Whole Blood	100/Rack
<b>BEC 366643*</b>	16 x 100		10.0 mL	Spray-Coated K2EDTA, Whole Blood	100/Rack
<b>BEC 366664</b>	13 x 75		2.0 mL	37 USP Units of Lithium Heparin, Plasma	100/Rack
<b>BEC 367884</b>	13 x 75		4.0 mL	60 USP Units of Lithium Heparin, Plasma	100/Rack
<b>BEC 367886</b>	13 x 100		6.0 mL	90 USP Units of Lithium Heparin, Plasma	100/Rack
<b>BEC 367871</b>	13 x 75		4.0 mL	60 USP Units of Sodium Heparin	100/Rack
<b>BEC 367878</b>	13 x 100		6.0 mL	90 USP Units of Sodium Heparin	100/Rack
<b>BEC 367874</b>	16 x 100		10.0 mL	150 USP Units of Sodium Heparin	100/Rack
<b>BEC 367921</b>	13 x 75		2.0 mL	Fluoride (Glucose), Plasma	100/Rack
<b>BEC 367922</b>	13 x 75		4.0 mL	Fluoride (Glucose), Plasma	100/Rack
<b>BEC 367925</b>	13 x 100		6.0 mL	Fluoride (Glucose), Plasma	100/Rack
<b>BEC 368056</b>	13 x 75		3.0 mL	PST™ Lithium Heparin w/ Gel Separator	100/Rack
<b>BEC 367962</b>	13 x 100		4.5 mL	PST™ Lithium Heparin w/ Gel Separator	100/Rack
<b>BEC 367964</b>	16 x 100		8.0 mL	PST™ Lithium Heparin w/ Gel Separator	100/Rack
<b>BEC 367981</b>	13 x 75		3.5 mL	SST™ Clot Activator/Polymer Gel, Silicon-Coated	100/Rack
<b>BEC 367983</b>	13 x 75		3.5 mL	SST™ Clot Activator/Polymer Gel, Silicon-Coated	100/Rack
<b>BEC 367977</b>	13 x 100		4.0 mL	SST™ Clot Activator/Polymer Gel, Silicon-Coated	100/Rack
<b>BEC 367986</b>	13 x 100		5.0 mL	SST™ Clot Activator/Polymer Gel, Silicon-Coated	100/Rack
<b>BEC 367988</b>	16 x 100		8.5 mL	SST™ Clot Activator/Polymer Gel, Silicon-Coated	100/Rack
<b>BEC 367985</b>	16 x 125		10.0 mL	SST™ Clot Activator/Double Polymer Gel, Silicon-Coated	100/Rack
<b>BEC 363083</b>	13 x 75		2.7 mL	Buffered NA Citrate (0.109M, 3.2%)	100/Rack
<b>BEC 368380</b>	13 x 100		6.0 mL	Serum Clot Activator, Silicon-Coated	100/Rack
<b>BEC 368381</b>	13 x 100		6.0 mL	K2EDTA	100/Rack
<b>BEC 367855</b>	13 x 75		3.0 mL	K2EDTA 5.4mg	100/Rack




**BD® ECLIPSE™ NEEDLE**

Bevel-oriented needle allows for low-angle injections

**BD® ECLIPSE™ NEEDLE W/ SMARTSLIP™ TECHNOLOGY**

Metal clip helps ensure a secure connection to luer slip syringes

Needle remains locked inside activated cover



**PRIOR TO ACTIVATION**

Wide-textured finger pad provides safe activation area for thumb and forefinger

**LUER LOCK SYRINGE**

Standard hub fits all Luer Lock syringes



**LUER SLIP SYRINGE**

**AFTER ACTIVATION**



**BD® ECLIPSE™ NEEDLES**

ITEM #	NEEDLE SIZE	QTY
<b>BEC 305761</b>	25G x 1"	100/Box
<b>BEC 305767</b>	25G x 1.5"	100/Box
<b>BEC 305762</b>	23G x 1"	100/Box
<b>BEC 305768</b>	22G x 1"	100/Box
<b>BEC 305763</b>	22G x 1.5"	100/Box
<b>BEC 305764</b>	21G x 1"	100/Box
<b>BEC 305765</b>	21G x 1.5"	100/Box
<b>BEC 305766</b>	18G x 1.5"	100/Box
<b>BEC 305759</b>	25G x 0.625"	100/Box
<b>BEC 305758</b>	27G x 0.5"	100/Box
<b>BEC 305757</b>	30G x 0.5"	100/Box

**BD® ECLIPSE™ NEEDLE WITH LUER-SLIP SYRINGE**

ITEM #	NEEDLE SIZE	QTY
<b>BEC 305770</b>	27G x 0.5"	100/Box
<b>BEC 305771</b>	30G x 0.5"	100/Box

**1ML LUER SLIP SYRINGE W/DETACHABLE NEEDLE**

ITEM #	NEEDLE SIZE	QTY
<b>BEC 305776</b>	27G x 0.5"	50/Box
<b>BEC 305777</b>	25G x 0.625"	50/Box

**BD® ECLIPSE™ NEEDLES WITH LUER-LOK™ SYRINGE**

**1ML BD® LUER-LOK™ SYRINGE W/DETACHABLE NEEDLE**

ITEM #	NEEDLE SIZE	QTY
<b>BEC 305778</b>	30G x 0.5"	300/Case
<b>BEC 305789</b>	27G x 0.5"	300/Case
<b>BEC 305780</b>	25G x 0.625"	300/Case

**3ML BD® LUER-LOK™ SYRINGE W/DETACHABLE NEEDLE**

ITEM #	NEEDLE SIZE	QTY
<b>BEC 305781</b>	25G x 0.625"	300/Case
<b>BEC 305787</b>	25G x 1"	300/Case
<b>BEC 305782</b>	23G x 1"	300/Case
<b>BEC 305788</b>	22G x 1"	300/Case
<b>BEC 305783</b>	22G x 1.5"	300/Case
<b>BEC 305779</b>	21G x 1"	300/Case
<b>BEC 305784</b>	21G x 1.5"	300/Case

**5ML BD® LUER-LOK™ SYRINGE W/DETACHABLE NEEDLE**

ITEM #	NEEDLE SIZE	QTY
<b>BEC 305785</b>	22G x 1.5"	300/Case

**10ML BD® LUER-LOK™ SYRINGE W/DETACHABLE NEEDLE**

ITEM #	NEEDLE SIZE	QTY
<b>BEC 305786</b>	22G x 1.5"	300/Case