



SINGLE URINE SPECIMEN COLLECTION KIT

item # TNR HXK44



- **Includes:** 90mL (3oz) Sterile Container w/Temperature Strip, Evidence Seal & Transport Bag w/Absorbent Pad
- Cup is made of polypropylene (cup) while lid is made of HDPE
- Packaged 100 Kits Per Case
- Kits are Individually Bagged



PHONE: 800.331.2716 | FAX: 800.359.8807 | WEB: MERCEDESMEDICAL.COM

All offers are subject to availability and are subject to change without notice. The names of various manufacturers, their instruments and their products referred to herein may be protected by trademark or other law, and are used herein solely for purpose of reference. Mercedes Medical expressly disclaims any affiliation with or sponsorship by the companies referred to herein. Tanner Scientific and all associated marks are trademarks of Tanner Scientific LLC. To stop receiving faxes, please fill out the information below and fax back to the number above. Name _____ Fax Number _____

67773
REV030918