



## SAFETY DATA SHEET In Accordance with ISO/DIS 11014

Review Date 02/20/2015

Revision Date 01/23/2013

### 1 Identification of the Substance/Mixture and Company

**Product Name:** Trans\*Form Complete

**MSDS Code:** SSCWTFR

**Product Description:** 10% Formalin Neutralizer

**Manufacturer/Supplier:**

American MasterTech  
1330 Thurman Street  
Lodi, CA 95240  
USA  
(800) 860-4073

**European Authorized Representative:**

Emergo Europe  
Moleensraat 15  
2513 BH The Hague,  
The Netherlands

**Emergency Telephone Number:** Infotrac (800) 535-5053 (24 hours) - International (011) 352-323-3500

### 2 Hazards Identification

**Classification of the substance or mixture**



GHS07

Skin Irrit. 2 H315 Causes skin irritation.

H319 Causes serious eye irritation.

**Classification according to Directive 67/548/EEC or Directive 1999/45/EC**



Xi; Irritant

R36/38: Irritating to eyes and skin.

**Information concerning particular hazards for human and environment:** Not applicable.

**Label elements**

**GHS label elements** The substance is classified and labeled according to the Globally Harmonized System (GHS).

**Hazard pictograms**



GHS07

**Signal word** Warning

**Hazard statements**

Causes skin irritation.

Causes serious eye irritation.

**Precautionary statements**

If medical advice is needed, have product container or label at hand.

Keep out of reach of children.

# SAFETY DATA SHEET

## In Accordance with ISO/DIS 11014

Review Date 02/20/2015

Revision Date 01/23/2013

**Product Name: Trans\*Form Complete**

Read label before use.

Wear protective gloves/protective clothing/eye protection/face protection.

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Specific treatment (see on this label).

Take off contaminated clothing and wash before reuse.

If skin irritation occurs: Get medical advice/attention.

If eye irritation persists: Get medical advice/attention.

**Classification system:****NFPA ratings (scale 0 - 4)**

Health = 1

Fire = 1

Reactivity = 1

**HMIS-ratings (scale 0 - 4)**

Health = 1

Fire = 1

Reactivity = 1

**Other hazards****Results of PBT and vPvB assessment**

**PBT:** Not applicable.

**vPvB:** Not applicable.

### 3 Composition

**Chemical characterization: Substances****CAS No. Description**

Proprietary Compound

**Identification number(s)**

### 4 First Aid Measures

**Description of first aid measures**

**After inhalation:** In case of unconsciousness place patient stably in side position for transportation.

**After skin contact:** Immediately wash with water and soap and rinse thoroughly.

**After eye contact:** Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.

**After swallowing:** If symptoms persist consult doctor.

**Information for doctor:**

**Most important symptoms and effects, both acute and delayed:** No further relevant information available.

**Indication of any immediate medical attention and special treatment needed:**

No further relevant information available.

### 5 Firefighting Measures

**Extinguishing media****Suitable extinguishing agents:**

CO<sub>2</sub>, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.

**Special hazards arising from the substance or mixture:** No further relevant information available.

# SAFETY DATA SHEET

## In Accordance with ISO/DIS 11014

Review Date 02/20/2015

Revision Date 01/23/2013

---

**Product Name: Trans\*Form Complete**


---

**Advice for firefighters****Protective equipment:** No special measures required.

## 6 Accidental Release Measures

**Personal precautions, protective equipment and emergency procedures:** Not required.**Environmental precautions:** No special measures required.**Methods and material for containment and cleaning up:** Pick up mechanically.**Reference to other sections**

See Section 7 for information on safe handling.

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

## 7 Handling and Storage

**Handling:****Precautions for safe handling:** No special precautions are necessary if used correctly.**Information about protection against explosions and fires:** No special measures required.**Conditions for safe storage, including any incompatibilities****Storage:****Requirements to be met by storerooms and receptacles:** No special requirements.**Information about storage in one common storage facility:** Not required.**Further information about storage conditions:** Keep receptacle tightly sealed.**Specific end use(s):** No further relevant information available.

## 8 Exposure Controls/Personal Protection

**Additional information about design of technical systems:** No further data; see item 7.**Control parameters****Components with limit values that require monitoring at the workplace:** Not required.**Additional information:** The lists that were valid during the creation were used as basis.**Exposure controls:****Personal protective equipment:****General protective and hygienic measures:**

Keep away from foodstuffs, beverages and feed.

Immediately remove all soiled and contaminated clothing.

Wash hands before breaks and at the end of work.

Avoid contact with the eyes and skin.

**Breathing equipment:** Not required.**Protection of hands:**

Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

# SAFETY DATA SHEET

## In Accordance with ISO/DIS 11014

Review Date 02/20/2015

Revision Date 01/23/2013

---

**Product Name: Trans\*Form Complete**


---

**Material of gloves:**

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer.

**Penetration time of glove material:**

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

**Eye protection:** Tightly sealed goggles

---

### 9 Physical and Chemical Properties

<b>Information on basic physical and chemical properties</b>	
<b>General Information</b>	
<b>Appearance</b>	
<b>Form:</b>	Solid
<b>Color:</b>	Not determined.
<b>Odor:</b>	Characteristic
<b>Odor threshold:</b>	Not determined.
<b>pH-value:</b>	Not applicable.
<b>Change in condition</b>	
<b>Melting point/Melting range:</b>	Undetermined.
<b>Boiling point/Boiling range:</b>	Undetermined.
<b>Flash point:</b>	Not applicable.
<b>Flammability (solid, gaseous):</b>	Product is not flammable.
<b>Ignition temperature:</b>	
<b>Decomposition temperature:</b>	Not determined.
<b>Auto igniting:</b>	Not determined.
<b>Danger of explosion:</b>	Product does not present an explosion hazard.
<b>Explosion limits:</b>	
<b>Lower:</b>	Not determined.
<b>Upper:</b>	Not determined.
<b>Vapor pressure:</b>	Not applicable.
<b>Density:</b>	Not determined.
<b>Relative density:</b>	Not determined.
<b>Vapour density:</b>	Not applicable.
<b>Evaporation rate:</b>	Not applicable.
<b>Solubility in / Miscibility with</b>	
<b>Water:</b>	Soluble.
<b>Segregation coefficient (n-octanol/water):</b> Not determined.	
<b>Dynamic:</b>	Not applicable.
<b>Kinematic:</b>	Not applicable.
<b>Other information:</b>	No further relevant information available.

# SAFETY DATA SHEET

## In Accordance with ISO/DIS 11014

Review Date 02/20/2015

Revision Date 01/23/2013

---

**Product Name: Trans\*Form Complete**


---

### 10 Stability and Reactivity

**Reactivity:****Chemical stability:****Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.**Possibility of hazardous reactions:** No dangerous reactions known.**Conditions to avoid:** No further relevant information available.**Incompatible materials:** No further relevant information available.**Hazardous decomposition products:** No dangerous decomposition products known.

### 11 Toxicology Information

**Information on toxicological effects****Acute toxicity:****Primary irritant effect:****On the skin:** Irritant to skin and mucous membranes.**On the eye:** Irritating effect.**Sensitization:** No sensitizing effects known.**Additional toxicological information:****Carcinogenic categories****IARC (International Agency for Research on Cancer)**

Substance is not listed.

**NTP (National Toxicology Program)**

Substance is not listed.

### 12 Ecological Information

**Toxicity:****Aquatic toxicity:** No further relevant information available.**Persistence and degradability:** No further relevant information available.**Behavior in environmental systems:****Bioaccumulative potential:** No further relevant information available.**Mobility in soil:** No further relevant information available.**Results of PBT and vPvB assessment****PBT:** Not applicable.**vPvB:** Not applicable.**Other adverse effects:** No further relevant information available.

### 13 Disposal Considerations

**Waste treatment methods:****Recommendation:** Disposal must be made according to official regulations.**Uncleaned packagings:****Recommendation:** Disposal must be made according to official regulations.

# SAFETY DATA SHEET

## In Accordance with ISO/DIS 11014

Review Date 02/20/2015

Revision Date 01/23/2013

**Product Name: Trans\*Form Complete**
**Recommended cleansing agent:** Water, if necessary with cleansing agents.

### 14 Transport Information

<b>UN-Number</b> DOT, ADR, ADN, IMDG, IATA	Void
<b>UN proper shipping name</b> DOT, ADR, ADN, IMDG, IATA	Void
<b>Transport hazard class(es)</b> DOT, ADR, ADN, IMDG, IATA Class	Void
<b>Packing group</b> DOT, ADR, IMDG, IATA	Void
<b>Environmental hazards</b> Marine pollutant:	No
<b>Special precautions for user:</b>	Not applicable.
<b>Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code:</b>	Not applicable.
<b>UN "Model Regulation":</b>	-

### 15 Regulatory Information

Safety, health and environmental regulations/legislation specific for the substance or mixture  
SARA

**Section 355 (extremely hazardous substances):**

Substance is not listed.

**Section 313 (Specific toxic chemical listings):**

Substance is not listed.

**TSCA (Toxic Substances Control Act):**

Substance is not listed.

**Proposition 65**
**Chemicals known to cause cancer:**

Substance is not listed.

**Chemicals known to cause reproductive toxicity for females:**

Substance is not listed.

**Chemicals known to cause reproductive toxicity for males:**

Substance is not listed.

**Chemicals known to cause developmental toxicity:**

Substance is not listed.

# SAFETY DATA SHEET

## In Accordance with ISO/DIS 11014

Review Date 02/20/2015

Revision Date 01/23/2013

---

**Product Name: Trans\*Form Complete**


---

**Carcinogenic categories****EPA (Environmental Protection Agency)**

Substance is not listed.

**TLV (Threshold Limit Value established by ACGIH)**

Substance is not listed.

**NIOSH-Ca (National Institute for Occupational Safety and Health)**

Substance is not listed.

**OSHA-Ca (Occupational Safety & Health Administration)**

Substance is not listed.

**GHS label elements** The substance is classified and labeled according to the Globally Harmonized System (GHS).**Hazard pictograms**

GHS07

**Signal word** Warning**Hazard statements**

Causes skin irritation.

Causes serious eye irritation.

**Precautionary statements**

If medical advice is needed, have product container or label at hand.

Keep out of reach of children.

Read label before use.

Wear protective gloves/protective clothing/eye protection/face protection.

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.

Continue rinsing.

Specific treatment (see on this label).

Take off contaminated clothing and wash before reuse.

If skin irritation occurs: Get medical advice/attention.

If eye irritation persists: Get medical advice/attention.

**Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

---

## 16 Other Information

This information contained herein is believed to be accurate but is not warranted to be so. Data and calculations are based on information furnished by the manufacturer of the product and manufacturers of the components of the product. Users are advised to confirm in advance of need that information is current, applicable and suited to the circumstances of use. Vendor assumes no responsibility for injury to vendee or third persons proximately caused by the material if reasonable safety procedures are not adhered to as stipulated in the data sheet. Furthermore, vendor assumes no responsibility for injury caused by abnormal use of the material even if reasonable safety procedures are followed. Any questions regarding this product should be directed to the manufacturer of the product as described in Section 1.

**Department issuing SDS:** Regulatory Department**Contact:** Phone (800) 860-4073

**SAFETY DATA SHEET**  
**In Accordance with ISO/DIS 11014**

Review Date 02/20/2015

Revision Date 01/23/2013

---

**Product Name: Trans\*Form Complete**

---

**Abbreviations and acronyms:**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

ACGIH: American Conference of Governmental Industrial Hygienists

CAS: Chemical Abstracts Service (division of the American Chemical Society)

NFPA: National Fire Protection Association (USA)

HMIS: Hazardous Materials Identification System (USA)

**Copyright © 2015 American MasterTech Scientific Incorporated. All rights reserved.**

No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.

---

— USA —