

03/10/2015

Kit Components

| Product code | Description |
|--------------|-------------|
|--------------|-------------|

| | |
|----------------|---|
| SSKT093 | Methenamine-Silver Compound, Sodium Borate |
|----------------|---|

Components:

| | |
|---------|-----------------------------|
| SSKC093 | Methenamine Silver Compound |
| SSKC095 | Sodium Borate Concentrate |



SAFETY DATA SHEET
In Accordance with ISO/DIS 11014

Review Date 03/10/2015

Revision Date 06/22/2012

1 Identification of the Substance/Mixture and Company

Product Name: Methenamine Silver Compound

MSDS Code: SSKC093

Product Description: Kit Component

Manufacturer/Supplier:

American MasterTech
1330 Thurman Street
Lodi, CA 95240
USA
(800) 860-4073

European Authorized Representative:

Emergo Europe
Moleensraat 15
2513 BH The Hague,
The Netherlands

Emergency Telephone Number: Infotrac (800) 535-5053 (24 hours) - International (011) 352-323-3500

2 Hazards Identification

Classification of the substance or mixture



GHS06 Skull and crossbones

Acute Tox. 3 H301 Toxic if swallowed.



GHS05 Corrosion

Skin Corr. 1B H314 Causes severe skin burns and eye damage.



GHS09 Environment

Aquatic Acute 1

Aquatic Chronic 1 H410 Very toxic to aquatic life with long lasting effects.



GHS07

Skin Sens. 1 H317 May cause an allergic skin reaction.

Classification according to Directive 67/548/EEC or Directive 1999/45/EC



C; Corrosive

R34: Causes burns.



Xi; Sensitizing

R43: May cause sensitisation by skin contact.

SAFETY DATA SHEET

In Accordance with ISO/DIS 11014

Review Date 03/10/2015

Revision Date 06/22/2012

Product Name: Methenamine Silver Compound



N; Dangerous for the environment

R50/53: Very toxic to aquatic organisms, may cause long-term adverse effects in the aquatic environment.

Information concerning particular hazards for human and environment:

The product has to be labelled due to the calculation procedure of international guidelines.

Classification system:

The classification was made according to the latest editions of international substances lists, and expanded upon from company and literature data.

Label elements**GHS label elements** The product is classified and labeled according to the Globally Harmonized System (GHS).**Hazard pictograms**

GHS05

GHS06

GHS09

Signal word Danger**Hazard-determining components of labeling:**

silver nitrate

methenamine

Hazard statements

Toxic if swallowed.

Causes severe skin burns and eye damage.

May cause an allergic skin reaction.

Very toxic to aquatic life with long lasting effects.

Precautionary statements

If medical advice is needed, have product container or label at hand.

Keep out of reach of children.

Read label before use.

IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician.

IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.

Continue rinsing.

Immediately call a POISON CENTER or doctor/physician.

Store locked up.

Dispose of contents/container in accordance with local/regional/national/international regulations.

Classification system:**NFPA ratings (scale 0 - 4)**

Health = 3

Fire = 0

Reactivity = 0

HMIS-ratings (scale 0 - 4)

Health = 3

Fire = 0

Reactivity = 0

Other hazards**Results of PBT and vPvB assessment****PBT:** Not applicable.

SAFETY DATA SHEET

In Accordance with ISO/DIS 11014

Review Date 03/10/2015

Revision Date 06/22/2012

Product Name: Methenamine Silver Compound
vPvB: Not applicable.

3 Composition

Chemical characterization: Mixtures
Description: Mixture of the substances listed below with nonhazardous additions.

| Dangerous components: | | |
|-------------------------------------|---|------|
| CAS: 7761-88-8 EINECS: 231-853-9 | silver nitrate ⚠ Ox. Sol. 2, H272; ⚠ Skin Corr. 1B, H314; ⚠ Aquatic Acute 1; Aquatic Chronic 1, H410 | <60% |
| CAS: 100-97-0 EINECS: 202-905-8 | methenamine ⚠ Flam. Liq. 1, H224; Flam. Sol. 2, H228; ⚠ Skin Sens. 1, H317 | <40% |

4 First Aid Measures

Description of first aid measures
General information: Immediately remove any clothing soiled by the product.

After inhalation:

Supply fresh air and to be sure call for a doctor.

In case of unconsciousness place patient stably in side position for transportation.

After skin contact: Immediately wash with water and soap and rinse thoroughly.

After eye contact: Rinse opened eye for several minutes under running water. Then consult a doctor.

After swallowing: Drink copious amounts of water and provide fresh air. Immediately call a doctor.

Information for doctor:
Most important symptoms and effects, both acute and delayed: No further relevant information available.

Indication of any immediate medical attention and special treatment needed:

No further relevant information available.

5 Firefighting Measures

Extinguishing media
Suitable extinguishing agents:

 CO₂, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.

Special hazards arising from the substance or mixture: No further relevant information available.

Advice for firefighters
Protective equipment: No special measures required.

6 Accidental Release Measures

Personal precautions, protective equipment and emergency procedures:

Wear protective equipment. Keep unprotected persons away.

Environmental precautions:

Do not allow product to reach sewage system or any water course.

Inform respective authorities in case of seepage into water course or sewage system.

Do not allow to enter sewers/ surface or ground water.

Methods and material for containment and cleaning up:

Use neutralizing agent.

SAFETY DATA SHEET

In Accordance with ISO/DIS 11014

Review Date 03/10/2015

Revision Date 06/22/2012

Product Name: Methenamine Silver Compound

Dispose contaminated material as waste according to item 13.
Ensure adequate ventilation.

Reference to other sections

See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.

7 Handling and Storage

Handling:**Precautions for safe handling:** Thorough dedusting.**Information about protection against explosions and fires:** No special measures required.**Conditions for safe storage, including any incompatibilities****Storage:****Requirements to be met by storerooms and receptacles:** No special requirements.**Information about storage in one common storage facility:** Not required.**Further information about storage conditions:** Keep receptacle tightly sealed.**Specific end use(s):** No further relevant information available.

8 Exposure Controls/Personal Protection

Additional information about design of technical systems: No further data; see item 7.**Control parameters****Components with limit values that require monitoring at the workplace:****7761-88-8 silver nitrate**

| | |
|-----|---------------------------------|
| PEL | 0.01 mg/m ³ as Ag |
| REL | 0.01 mg/m ³ as Ag |
| TLV | 0.01 mg/m ³ as Ag |

Additional information: The lists that were valid during the creation were used as basis.**Exposure controls:****Personal protective equipment:****General protective and hygienic measures:**

Keep away from foodstuffs, beverages and feed.
Immediately remove all soiled and contaminated clothing.
Wash hands before breaks and at the end of work.
Avoid contact with the eyes and skin.

Breathing equipment: Not required.**Protection of hands:**

Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

SAFETY DATA SHEET

In Accordance with ISO/DIS 11014

Review Date 03/10/2015

Revision Date 06/22/2012

Product Name: Methenamine Silver Compound

Material of gloves:

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

Penetration time of glove material:

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

Eye protection: Tightly sealed goggles

9 Physical and Chemical Properties

Information on basic physical and chemical properties
General Information
Appearance

| | |
|------------------------|-----------------|
| Form: | Solid |
| Color: | Brown |
| Odor: | Nearly odorless |
| Odor threshold: | Not determined. |

pH-value: Not applicable.

Change in condition

| | |
|-------------------------------------|---------------|
| Melting point/Melting range: | Undetermined. |
| Boiling point/Boiling range: | Undetermined. |

Flash point: Not applicable.

Flammability (solid, gaseous): Not determined.

Ignition temperature:

Decomposition temperature: Not determined.

Auto igniting: Product is not selfigniting.

Danger of explosion: Product does not present an explosion hazard.

Explosion limits:

| | |
|---------------|-----------------|
| Lower: | Not determined. |
| Upper: | Not determined. |

Vapor pressure: Not applicable.

| | |
|--------------------------|-----------------|
| Density: | Not determined. |
| Relative density: | Not determined. |
| Vapour density: | Not applicable. |
| Evaporation rate: | Not applicable. |

Solubility in / Miscibility with

Water: Soluble.

Segregation coefficient (n-octanol/water): Not determined.

Dynamic: Not applicable.

Kinematic: Not applicable.

Other information: No further relevant information available.

SAFETY DATA SHEET

In Accordance with ISO/DIS 11014

Review Date 03/10/2015

Revision Date 06/22/2012

Product Name: Methenamine Silver Compound

10 Stability and Reactivity

Reactivity:**Chemical stability:****Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.**Possibility of hazardous reactions:** No dangerous reactions known.**Conditions to avoid:** No further relevant information available.**Incompatible materials:** No further relevant information available.**Hazardous decomposition products:** No dangerous decomposition products known.

11 Toxicology Information

Information on toxicological effects**Acute toxicity:****LD/LC50 values that are relevant for classification:****7761-88-8 silver nitrate**

Oral | LD50 | 50 mg/kg (mouse)

Primary irritant effect:**On the skin:** Caustic effect on skin and mucous membranes.**On the eye:** Strong caustic effect.**Sensitization:** Sensitization possible through skin contact.**Additional toxicological information:**

The product shows the following dangers according to internally approved calculation methods for preparations:

Corrosive

Irritant

Swallowing will lead to a strong caustic effect on mouth and throat and to the danger of perforation of esophagus and stomach.

Carcinogenic categories**IARC (International Agency for Research on Cancer)**

None of the ingredients are listed.

NTP (National Toxicology Program)

None of the ingredients are listed.

12 Ecological Information

Toxicity:**Aquatic toxicity:** No further relevant information available.**Persistence and degradability:** No further relevant information available.**Behavior in environmental systems:****Bioaccumulative potential:** No further relevant information available.**Mobility in soil:** No further relevant information available.**Ecotoxicological effects:****Remark:** Very toxic for fish**Additional ecological information:****General notes:**

Water hazard class 3 (Self-assessment): extremely hazardous for water

Do not allow product to reach ground water, water course or sewage system, even in small quantities.

SAFETY DATA SHEET

In Accordance with ISO/DIS 11014

Review Date 03/10/2015

Revision Date 06/22/2012

Product Name: Methenamine Silver Compound

Must not reach bodies of water or drainage ditch undiluted or unneutralized.
 Danger to drinking water if even extremely small quantities leak into the ground.
 Also poisonous for fish and plankton in water bodies.
 Very toxic for aquatic organisms
Results of PBT and vPvB assessment
PBT: Not applicable.
vPvB: Not applicable.
Other adverse effects: No further relevant information available.

13 Disposal Considerations

Uncleaned packagings:
Recommendation: Disposal must be made according to official regulations.
Recommended cleansing agent: Water, if necessary with cleansing agents.

14 Transport Information

| | |
|---|-----------------|
| UN-Number DOT, ADR, ADN, IMDG, IATA | Void |
| UN proper shipping name DOT, ADR, ADN, IMDG, IATA | Void |
| Transport hazard class(es) DOT, ADR, ADN, IMDG, IATA Class | Void |
| Packing group DOT, ADR, IMDG, IATA | Void |
| Environmental hazards Marine pollutant: | No |
| Special precautions for user: | Not applicable. |
| Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code: | Not applicable. |
| UN "Model Regulation": | - |

15 Regulatory Information

Safety, health and environmental regulations/legislation specific for the substance or mixture
SARA

Section 355 (extremely hazardous substances):

None of the ingredients are listed.

Section 313 (Specific toxic chemical listings):

| | |
|-----------|----------------|
| 7761-88-8 | silver nitrate |
|-----------|----------------|

SAFETY DATA SHEET

In Accordance with ISO/DIS 11014

Review Date 03/10/2015

Revision Date 06/22/2012

Product Name: Methenamine Silver Compound

| |
|--|
| TSCA (Toxic Substances Control Act): |
| All ingredients are listed. |
| Proposition 65 |
| Chemicals known to cause cancer: |
| None of the ingredients is listed. |
| Chemicals known to cause reproductive toxicity for females: |
| None of the ingredients are listed. |
| Chemicals known to cause reproductive toxicity for males: |
| None of the ingredients are listed. |
| Chemicals known to cause developmental toxicity: |
| None of the ingredients are listed. |

Carcinogenic categories

| |
|---|
| EPA (Environmental Protection Agency) |
| None of the ingredients are listed. |
| TLV (Threshold Limit Value established by ACGIH) |
| None of the ingredients are listed. |
| NIOSH-Ca (National Institute for Occupational Safety and Health) |
| None of the ingredients are listed. |
| OSHA-Ca (Occupational Safety & Health Administration) |
| None of the ingredients are listed. |

GHS label elements The product is classified and labeled according to the Globally Harmonized System (GHS).

Hazard pictograms


GHS05 GHS06 GHS09

Signal word Danger

Hazard-determining components of labelling:

 silver nitrate
 methenamine

Hazard statements

 Toxic if swallowed.
 Causes severe skin burns and eye damage.
 May cause an allergic skin reaction.
 Very toxic to aquatic life with long lasting effects.

Precautionary statements

 If medical advice is needed, have product container or label at hand.
 Keep out of reach of children.
 Read label before use.

IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician.

IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.

 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.
 Continue rinsing.

SAFETY DATA SHEET

In Accordance with ISO/DIS 11014

Review Date 03/10/2015

Revision Date 06/22/2012

Product Name: Methenamine Silver Compound

Immediately call a POISON CENTER or doctor/physician.

Store locked up.

Dispose of contents/container in accordance with local/regional/national/international regulations.

Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

16 Other Information

This information contained herein is believed to be accurate but is not warranted to be so. Data and calculations are based on information furnished by the manufacturer of the product and manufacturers of the components of the product. Users are advised to confirm in advance of need that information is current, applicable and suited to the circumstances of use. Vendor assumes no responsibility for injury to vendee or third persons proximately caused by the material if reasonable safety procedures are not adhered to as stipulated in the data sheet. Furthermore, vendor assumes no responsibility for injury caused by abnormal use of the material even if reasonable safety procedures are followed. Any questions regarding this product should be directed to the manufacturer of the product as described in Section 1.

Department issuing SDS: Regulatory Department

Contact: Phone (800) 860-4073

Abbreviations and acronyms:

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

ACGIH: American Conference of Governmental Industrial Hygienists

NFPA: National Fire Protection Association (USA)

HMIS: Hazardous Materials Identification System (USA)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

Copyright © 2015 American MasterTech Scientific Incorporated. All rights reserved.

No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.



SAFETY DATA SHEET In Accordance with ISO/DIS 11014

Review Date 03/10/2015

Revision Date 06/21/2012

1 Identification of the Substance/Mixture and Company

Product Name: Sodium Borate Concentrate

MSDS Code: SSKC095

Product Description: Kit Component

Manufacturer/Supplier:

American MasterTech
1330 Thurman Street
Lodi, CA 95240
USA
(800) 860-4073

European Authorized Representative:

Emergo Europe
Moleensraat 15
2513 BH The Hague,
The Netherlands

Emergency Telephone Number: Infotrac (800) 535-5053 (24 hours) - International (011) 352-323-3500

2 Hazards Identification

Classification of the substance or mixture



GHS08 Health hazard

Repr. 1B H360 May damage fertility or the unborn child.

Classification according to Directive 67/548/EEC or Directive 1999/45/EC Not applicable.

Information concerning particular hazards for human and environment:

The product has to be labelled due to the calculation procedure of international guidelines.

Classification system:

The classification was made according to the latest editions of international substances lists, and expanded upon from company and literature data.

Label elements

GHS label elements The product is classified and labeled according to the Globally Harmonized System (GHS).

Hazard pictograms



GHS08

Signal word Danger

Hazard statements

May damage fertility or the unborn child.

Precautionary statements

If medical advice is needed, have product container or label at hand.

Keep out of reach of children.

Read label before use.

Use personal protective equipment as required.

SAFETY DATA SHEET

In Accordance with ISO/DIS 11014

Review Date 03/10/2015

Revision Date 06/21/2012

Product Name: Sodium Borate Concentrate

Obtain special instructions before use.

Do not handle until all safety precautions have been read and understood.

IF exposed or concerned: Get medical advice/attention.

Store locked up.

Dispose of contents/container in accordance with local/regional/national/international regulations.

Classification system:

NFPA ratings (scale 0 - 4)

Health = 1

Fire = 0

Reactivity = 0

HMIS-ratings (scale 0 - 4)

Health = 1

Fire = 0

Reactivity = 0

Other hazards

Results of PBT and vPvB assessment

PBT: Not applicable.


vPvB: Not applicable.

3 Composition

Chemical characterization: Mixtures

Description: Mixture of the substances listed below with nonhazardous additions.

Dangerous components:

| | | | |
|-------------------|-----------------------------------|--|------|
| CAS: 1303-96-4 | Disodium tetraborate, decahydrate |  Repr. 1B, H360 | <10% |
| EINECS: 215-540-4 | | | |

SVHC

| | |
|-----------|-----------------------------------|
| 1303-96-4 | Disodium tetraborate, decahydrate |
|-----------|-----------------------------------|

4 First Aid Measures

Description of first aid measures

After inhalation: Supply fresh air; consult doctor in case of complaints.

After skin contact: Generally the product does not irritate the skin.

After eye contact: Rinse opened eye for several minutes under running water.

After swallowing: If symptoms persist consult doctor.

Information for doctor:

Most important symptoms and effects, both acute and delayed: No further relevant information available.

Indication of any immediate medical attention and special treatment needed:

No further relevant information available.

5 Firefighting Measures

Extinguishing media

Suitable extinguishing agents:

CO₂, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.

Special hazards arising from the substance or mixture: No further relevant information available.

SAFETY DATA SHEET

In Accordance with ISO/DIS 11014

Review Date 03/10/2015

Revision Date 06/21/2012

Product Name: Sodium Borate Concentrate

Advice for firefighters**Protective equipment:** No special measures required.

6 Accidental Release Measures

Personal precautions, protective equipment and emergency procedures: Not required.**Environmental precautions:**

Dilute with plenty of water.

Do not allow to enter sewers/ surface or ground water.

Methods and material for containment and cleaning up:

Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).

Reference to other sections

See Section 7 for information on safe handling.

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

7 Handling and Storage

Handling:**Precautions for safe handling:** No special precautions are necessary if used correctly.**Information about protection against explosions and fires:** No special measures required.**Conditions for safe storage, including any incompatibilities****Storage:****Requirements to be met by storerooms and receptacles:** No special requirements.**Information about storage in one common storage facility:** Not required.**Further information about storage conditions:** None.**Specific end use(s):** No further relevant information available.

8 Exposure Controls/Personal Protection

Additional information about design of technical systems: No further data; see item 7.**Control parameters****Components with limit values that require monitoring at the workplace:****1303-96-4 Disodium tetraborate, decahydrate**

| | |
|-----|---|
| REL | 5 mg/m ³ |
| TLV | Short-term value: 6* mg/m ³ Long-term value: 2* mg/m ³ *as inhalable fraction |

Additional information: The lists that were valid during the creation were used as basis.**Exposure controls:****Personal protective equipment:****General protective and hygienic measures:** Wash hands before breaks and at the end of work.**Breathing equipment:** Not required.**Protection of hands:**

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

SAFETY DATA SHEET

In Accordance with ISO/DIS 11014

Review Date 03/10/2015

Revision Date 06/21/2012

Product Name: Sodium Borate Concentrate

Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

Material of gloves:

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

Penetration time of glove material:

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

Eye protection: Goggles recommended during refilling.

9 Physical and Chemical Properties

| | |
|--|---|
| Information on basic physical and chemical properties | |
| General Information | |
| Appearance | |
| Form: | Fluid |
| Color: | Colorless |
| Odor: | Odorless |
| Odor threshold: | Not determined. |
| pH-value: | Not determined. |
| Change in condition | |
| Melting point/Melting range: | Undetermined. |
| Boiling point/Boiling range: | 100°C (212 °F) |
| Flash point: | Not applicable. |
| Flammability (solid, gaseous): | Not applicable. |
| Ignition temperature: | |
| Decomposition temperature: | Not determined. |
| Auto igniting: | Product is not selfigniting. |
| Danger of explosion: | Product does not present an explosion hazard. |
| Explosion limits: | |
| Lower: | Not determined. |
| Upper: | Not determined. |
| Vapor pressure at 20°C (68 °F): | 23 hPa (17 mm Hg) |
| Density: | Not determined. |
| Relative density: | Not determined. |
| Vapour density: | Not determined. |
| Evaporation rate: | Not determined. |
| Solubility in / Miscibility with | |
| Water: | Fully miscible. |
| Segregation coefficient (n-octanol/water): | |
| Dynamic: | Not determined. |
| Kinematic: | Not determined. |

SAFETY DATA SHEET

In Accordance with ISO/DIS 11014

Review Date 03/10/2015

Revision Date 06/21/2012

Product Name: Sodium Borate Concentrate

| | |
|---------------------------|--|
| Other information: | No further relevant information available. |
|---------------------------|--|

10 Stability and Reactivity

Reactivity:**Chemical stability:****Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.**Possibility of hazardous reactions:** No dangerous reactions known.**Conditions to avoid:** No further relevant information available.**Incompatible materials:** No further relevant information available.**Hazardous decomposition products:** No dangerous decomposition products known.

11 Toxicology Information

Information on toxicological effects**Acute toxicity:****Primary irritant effect:****On the skin:** No irritant effect.**On the eye:** No irritating effect.**Sensitization:** No sensitizing effects known.**Additional toxicological information:****Carcinogenic categories****IARC (International Agency for Research on Cancer)**

None of the ingredients are listed.

NTP (National Toxicology Program)

None of the ingredients are listed.

12 Ecological Information

Toxicity:**Aquatic toxicity:** No further relevant information available.**Persistence and degradability:** No further relevant information available.**Behavior in environmental systems:****Bioaccumulative potential:** No further relevant information available.**Mobility in soil:** No further relevant information available.**Additional ecological information:****General notes:**

Water hazard class 1 (Self-assessment): slightly hazardous for water

Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

Results of PBT and vPvB assessment**PBT:** Not applicable.**vPvB:** Not applicable.**Other adverse effects:** No further relevant information available.

SAFETY DATA SHEET

In Accordance with ISO/DIS 11014

Review Date 03/10/2015

Revision Date 06/21/2012

Product Name: Sodium Borate Concentrate

13 Disposal Considerations

Uncleaned packagings:**Recommendation:** Disposal must be made according to official regulations.**Recommended cleansing agent:** Water, if necessary with cleansing agents.

14 Transport Information

| | |
|---|-----------------|
| UN-Number DOT, ADR, ADN, IMDG, IATA | Void |
| UN proper shipping name DOT, ADR, ADN, IMDG, IATA | Void |
| Transport hazard class(es) DOT, ADR, ADN, IMDG, IATA Class | Void |
| Packing group DOT, ADR, IMDG, IATA | Void |
| Environmental hazards Marine pollutant: | No |
| Special precautions for user: | Not applicable. |
| Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code: | Not applicable. |
| UN "Model Regulation": | - |

15 Regulatory Information

Safety, health and environmental regulations/legislation specific for the substance or mixture
SARA

Section 355 (extremely hazardous substances):

None of the ingredients are listed.

Section 313 (Specific toxic chemical listings):

None of the ingredients is listed.

TSCA (Toxic Substances Control Act):

All ingredients are listed.

Proposition 65**Chemicals known to cause cancer:**

None of the ingredients is listed.

Chemicals known to cause reproductive toxicity for females:

None of the ingredients are listed.

Chemicals known to cause reproductive toxicity for males:

None of the ingredients are listed.

SAFETY DATA SHEET

In Accordance with ISO/DIS 11014

Review Date 03/10/2015

Revision Date 06/21/2012

Product Name: Sodium Borate Concentrate

| | |
|---|-----------------------------------|
| Chemicals known to cause developmental toxicity: | |
| None of the ingredients are listed. | |
| Carcinogenic categories | |
| EPA (Environmental Protection Agency) | |
| None of the ingredients are listed. | |
| TLV (Threshold Limit Value established by ACGIH) | |
| 1303-96-4 | Disodium tetraborate, decahydrate |
| | A4 |
| NIOSH-Ca (National Institute for Occupational Safety and Health) | |
| None of the ingredients are listed. | |
| OSHA-Ca (Occupational Safety & Health Administration) | |
| None of the ingredients are listed. | |

GHS label elements The product is classified and labeled according to the Globally Harmonized System (GHS).

Hazard pictograms


GHS08

Signal word Danger

Hazard statements

May damage fertility or the unborn child.

Precautionary statements

If medical advice is needed, have product container or label at hand.

Keep out of reach of children.

Read label before use.

Use personal protective equipment as required.

Obtain special instructions before use.

Do not handle until all safety precautions have been read and understood.

IF exposed or concerned: Get medical advice/attention.

Store locked up.

Dispose of contents/container in accordance with local/regional/national/international regulations.

Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

16 Other Information

This information contained herein is believed to be accurate but is not warranted to be so. Data and calculations are based on information furnished by the manufacturer of the product and manufacturers of the components of the product. Users are advised to confirm in advance of need that information is current, applicable and suited to the circumstances of use. Vendor assumes no responsibility for injury to vendee or third persons proximately caused by the material if reasonable safety procedures are not adhered to as stipulated in the data sheet. Furthermore, vendor assumes no responsibility for injury caused by abnormal use of the material even if reasonable safety procedures are followed. Any questions regarding this product should be directed to the manufacturer of the product as described in Section 1.

Department issuing SDS: Regulatory Department

Contact: Phone (800) 860-4073

SAFETY DATA SHEET
In Accordance with ISO/DIS 11014

Review Date 03/10/2015

Revision Date 06/21/2012

Product Name: Sodium Borate Concentrate

Abbreviations and acronyms:

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

ACGIH: American Conference of Governmental Industrial Hygienists

NFPA: National Fire Protection Association (USA)

HMIS: Hazardous Materials Identification System (USA)

Copyright © 2015 American MasterTech Scientific Incorporated. All rights reserved.

No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.

— USA —