



SAFETY DATA SHEET

Doc. ID: A16792-75 Rev. AG
Revised (year/month/day) 2015/04/14

Section 1 Identification of the Substance/mixture and of the Company/undertaking

1.1 Product Identifier

Product Name Wash Buffer II
Part Number A16792
Series Name ACCESS®

1.2 Relevant identified uses of the substance or mixture and uses advised against

Product Use For In Vitro Diagnostic Use. See product literature for details.

1.3 Details of the supplier of the safety data sheet

Manufacturer

Beckman Coulter, Inc.
250 S. Kraemer Blvd
Brea, CA 92821, U.S.A.
Tel: 800-854-3633

EC REP Address

Beckman Coulter Eurocenter S.A.
22, rue Juste-Oliver, Case Postale 1044,
CH-1260 Nyon 1, Switzerland.
Telephone +41 (0)22 365 36 11
Monday through Friday, 9:00 am to
7:00pm)

e-mail address SDSNT@beckman.com

1.4 Emergency telephone number

Telephone number (24H) Chemtrec Emergency Tel No. U.S.A. 800-424-9300, International (001) 703-527-3887

Distributor and Emergency Phone No.

Refer to attached list, Document ID: [472050](#), for local distributor and emergency phone numbers.

Section 2 Hazards Identification

2.1 Classification of substance or mixture

Product Description Mixture
Colorless to pale yellow; Clear; Liquid; Odorless

Classification according to EC 1272/2008 (CLP/GHS)
Skin Sensitization Category 1

Classification according to EC Directives 1999/45/EC and 67/548/EEC
Xi;R43

Classification according to US-OSHA (HCS 29 CFR 1910.1200) and UN GHS
Not classified as hazardous per US-OSHA HCS 2012 and UN GHS

Section 2 Hazards Identification (Continued)

2.2 Label Elements

According to EC 1272/2008 (CLP/GHS), US-OSHA and UN GHS

Hazardous Ingredients

reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC# 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC# 220-239-6](3:1)

Pictogram



Signal Word

WARNING

Hazard Statements

H317 May cause an allergic skin reaction.

Precautionary Statements

- P261 Avoid breathing vapours.
 - P272 Contaminated work clothing should not be allowed out of the workplace.
 - P280 Wear protective gloves, protective clothing and eye/face protection.
 - P302+P352 IF ON SKIN: Wash with plenty of soap and water.
 - P333+P313 If skin irritation or rash occurs: Get medical advice/attention.
 - P362+P364 Take off contaminated clothing and wash it before use.
 - P501 Dispose of contents/container in accordance with local/national regulations
- Product label will display most significant precautionary statements.

2.3 Other hazards

Results of PBT and vPvB assessment

PBT: Not applicable.

vPvB: Not applicable.

This product contains concentrations of azide below the hazardous level which with repeated contact with lead and copper commonly found in plumbing drains may result in the build up of shock sensitive compounds. Sodium azide forms explosive compounds with heavy metals.

See Section 11 Toxicological Information for more detailed health information.

Section 3 Composition and Information on Ingredients

3.2 Mixtures

Hazardous Ingredients:		Hazard Classification of Pure Ingredients			
Chemical Name	% by wt.	EU-67/548/EEC	EU 1272/2008 CLP/GHS	GHS	
Sodium Azide CAS # 26628-22-8 EINECS # 247-852-1 Index # 011-004-00-7	< 0.1	T+;R28-32 N;R50/53	Acute Tox. Oral 2 Aquatic Acute 1 Aquatic Longterm 1 H300; H400; H410	Acute Tox. Oral 2 Aquatic Acute 1 Aquatic Longterm 1 H300; H400; H410	2, 8

SAFETY DATA SHEET

Doc. ID: A16792-75 AG
Revised (year/month/day) 2015/04/14

Section 3 Composition and Information on Ingredients (Continued)

reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC# 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC# 220-239-6](3:1) CAS # 55965-84-9 EINECS # Not available Index # Not available	< 0.05	T;R23/24/25-34-43 N;R50/53	Acute Tox. Dermal 3 Acute Tox. Inhal. 3 Acute Tox. Oral 3 Aquatic Acute 1 Aquatic Longterm 1 Skin Corr. 1B Skin Sens. 1 H301; H311; H314; H317; H331; H400; H410	Acute Tox. Dermal 3 Acute Tox. Inhal. 3 Acute Tox. Oral 3 Aquatic Acute 1 Aquatic Longterm 1 Skin Corr. 1B Skin Sens. 1 H301; H311; H314; H317; H331; H400; H410	9
--	--------	-------------------------------	---	---	---

2 - Substance with Community workplace exposure limits

8 - Present at concentration below the cut-off limits.

9 - Mixture of 5-chloro-2-methyl-4-isothiazolin-3-one [EC# 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC# 220-239-6] (3:1) is the active ingredient of ProClin 300.

See section 8 for available Occupational exposure limits

See Section 15 for additional regulatory information

See Section 16 for hazard class, hazard statements and risk phrase description

Section 4 First Aid Measures

4.1 Description of first aid measures

Inhalation

If product is inhaled, move exposed individual to fresh air. If individual is not breathing, begin artificial respiration immediately and obtain medical attention.

Eye Contact

If product enters eyes, wash eyes gently under running water for 15 minutes or longer, making sure that the eyelids are held open. If pain or irritation occur, obtain medical attention.

Skin Contact

In case of skin contact, flush with copious amounts of water for at least 15 minutes. If pain or irritation occur, obtain medical attention.

Ingestion

If ingested, wash mouth out with water. If irritation or discomfort occurs, seek medical attention.

4.2 Most important symptoms and effects, both acute and delayed

May produce an allergic reaction in some people.

See Section 11 Toxicological Information for more detailed health information.

4.3 Indication of any immediate medical attention and special treatment needed

No specific medical attention or treatment required.

Section 5 Fire Fighting Measures

Flammable Properties

Nonflammable aqueous solution.

5.1 Extinguishing Media

In case of fire use carbon dioxide (CO₂), dry chemical, water spray or foam. For large fires use extinguishing media suitable for surrounding fire.

5.2 Special hazards arising from the substance or mixture

Special Fire and Explosion Hazards

No special hazards determined.

Section 5 Fire Fighting Measures (Continued)

Hazardous Combustion Products

No combustion products posing significant hazards are expected from this product (an aqueous solution).

5.3 Advice for fire fighters

Protective Equipment

Self-contained breathing apparatus is recommended for firefighters in all chemical fire situations.

5.4 Additional information

No further relevant information available.

Section 6 Accidental Release Measures

6.1 Personal precautions, protective equipment and emergency procedures

Personal Precautions

Observe general safety guidelines for protection; avoid eye and skin contact. Wear protective gloves, protective clothing and eye/face protection.

6.2 Environmental Precautions

Contain spill to prevent migration.
Do not allow the undiluted product to enter sewers/surface or ground water.
Dispose of contents/container in accordance with local regulations

6.3 Methods and material for containment and cleaning up

Spill and Leak Procedures

Absorb spilled material with an appropriate inert, non-flammable absorbent and dispose according to local regulations.

6.4 Reference to other sections

Refer sections 8 and 13.

Section 7 Handling and Storage

7.1 Precautions for safe handling

Use good laboratory procedures; avoid eye and skin contact.

7.2 Conditions for safe storage, including any incompatibilities

Store at 15 to 30°C, as directed on the product label.
To maintain product quality, store according to the instructions in the product labeling.
Store away from strong acids, strong bases, strong oxidizers and incompatible materials (section 10).

7.3 Specific end uses

No further relevant information available.

Section 8 Exposure Controls and Personal Protection

8.1 Control parameters

Exposure Limits

US OSHA

None established

ACGIH

Sodium Azide
CAS # 26628-22-8

0.29 mg/m3 Ceiling (as NaN₃); 0.11 ppm Ceiling (as Hydrazoic acid) (vapor)

Section 8 Exposure Controls and Personal Protection (Continued)

DFG MAK

Sodium Azide
CAS # 26628-22-8

0.4 mg/m³ Peak (inhalable fraction); 0.2 mg/m³ TWA MAK (inhalable fraction)

Ireland

Sodium Azide
CAS # 26628-22-8

0.1 mg/m³ TWA (as NaN₃); 0.3 mg/m³ STEL (as NaN₃); Potential for cutaneous absorption

IOELVs

Sodium Azide
CAS # 26628-22-8

Possibility of significant uptake through the skin; 0.1 mg/m³ TWA; 0.3 mg/m³ STEL

NIOSH

None established

Japan

None established

8.2 Exposure controls

Engineering Controls

No special engineering controls are required. Use with good general ventilation.

Eye Protection

Safety glasses or chemical goggles should be worn to prevent eye contact. Refer U.S. OSHA 29 CFR 1910.133, European Standard EN166 or appropriate government standards.

Skin Protection

Impervious gloves, such as Nitrile or equivalent, should be worn to prevent skin contact. Refer U.S. OSHA 29 CFR 1910.138, European Standard EN374 or appropriate government standards.

Respiratory Protection

Under normal conditions, the use of this product should not require respiratory protection. If overexposure should occur and ventilation is not adequate to maintain airborne concentrations at acceptable levels, the use of respiratory protection should be evaluated by a qualified professional.

Section 9 Physical and Chemical Properties

9.1 Information on basic physical and chemical properties

Physical State	Liquid	Specific Gravity (Water=1.0)	≈ 1.01 @20°C
Color	Colorless to pale yellow	Solubility	
Transparency	Clear	Water	Miscible
Odor	Odorless	Organic	Not determined
pH	8.3	Partition coefficient: n-octanol/water	Not determined
Freezing Point	Not determined	Auto-ignition Temp.	Not applicable
Boiling Point	Not determined	Decomposition Temperature	Not determined
Flash Point	Not applicable	Percent Volatiles	Not applicable

Section 9 Physical and Chemical Properties (Continued)

Evaporation Rate	Not determined	Vapor Pressure	Not determined
Flammability (Solid, Gas)	Not applicable	Viscosity	Not determined
Flammability Limits	Not applicable	Explosive Properties	Not applicable
Vapor Density	Not determined	Oxidizing Properties	Not applicable
Odor Threshold	Not applicable		
9.2 Other Information	No further relevant information available.		

Section 10 Stability and Reactivity

10.1 Reactivity	No further relevant information available.
10.2 Chemical Stability	The product is stable in accordance with recommended storage conditions.
10.3 Possibility of hazardous reactions	Sodium azide forms explosive compounds with heavy metals. Repeated contact of low concentrations of azide with lead and copper commonly found in plumbing drains may result in the build up of shock sensitive compounds.
10.4 Conditions to Avoid	Avoid contact with incompatible materials. Avoid exposure to heat and direct sunlight.
10.5 Incompatible materials	Metals and metallic compounds
10.6 Hazardous Decomposition Products	No decomposition products posing significant hazards would be expected from this product (an aqueous solution).

Section 11 Toxicological Information

11.1 Information on toxicological effects	
Toxicity Data for Hazardous Ingredients	
Sodium Azide CAS # 26628-22-8	Oral LD50 Rat 27 mg/kg; Dermal LD50 Rat 50 mg/kg; Dermal LD50 Rabbit 20 mg/kg
Primary Routes of Exposure	Eye contact, ingestion, inhalation, and skin contact.
Skin Corrosion/Irritation	No data available.
Serious eye damage/eye irritation	No data available.
Respiratory/skin sensitization	May cause sensitization by skin contact.
Carcinogenicity	This product does not contain a reportable concentration ($\geq 0.1\%$) of any ingredient listed as carcinogen by ACGIH, IARC, NTP, OSHA or 1272/2008 EC regulation.
Germ cell mutagenicity	No data available.
Reproductive Toxicity	No data available.

Section 11 Toxicological Information (Continued)

Specific target organ toxicity – single exposure

No data available.

Specific target organ toxicity – repeated exposure

No data available.

Aspiration hazard

No data available.

Other Information

No further relevant information available.

Section 12 Ecological Information

12.1 Ecotoxicity

Fresh Water Species

Sodium Azide
CAS # 26628-22-8

96 h LC50 *Oncorhynchus mykiss*: 0.8 mg/L; 96 h LC50 *Lepomis macrochirus*:
0.7 mg/L; 96 h LC50 *Pimephales promelas*: 5.46 mg/L [flow-through]

Microtox

No information available.

Water Flea

No information available.

Fresh Water Algae

No information available.

12.2 Persistence and degradability Not determined for the product.

12.3 Bioaccumulation Not determined for the product.

12.4 Mobility in soil Not determined for the product.

12.5 Results of PBT and vPvB assessment

Not determined for the product. PBT: Not applicable, vPvB: Not applicable.

12.6 Other Adverse Effects

This product contains environmentally hazardous substance below the cutoff level. Refer section 3 for ingredient information. Do not allow undiluted product to enter sewer/surface or ground water.

Section 13 Disposal Considerations

13.1 Waste treatment methods

Product Waste Disposal

Chemical residues and remains should be routinely handled as special waste. This must be disposed of in compliance with anti-pollution and other laws of the country concerned. To ensure compliance we recommend that you contact the relevant (local) authorities and/or an approved waste-disposal company for information.

Sodium azide preservative may form explosive compounds in metal drain lines. See NIOSH Bulletin: Explosive Azide Hazard (8/16/76).

To avoid the possible build-up of azide compounds, flush wastepipes with water after the disposal of undiluted reagent. Sodium azide disposal must be in accordance with appropriate local regulations.

Package disposal

Dispose of waste product, unused product and contaminated packaging in compliance with federal, state and local regulations. If unsure of the applicable requirements, contact the authorities for information.

13.2 Additional information

Suggested European waste catalogue 18 01 07 - chemicals other than those mentioned in 18 01 06. Dispose in accordance with national, state and local waste regulations.

Section 14 Transport Information

Transportation of this product is not regulated under ICAO, IMDG, US DOT, European ADR or Canadian TDG.

Section 15 Regulatory Information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

US Federal and State Regulations

SARA 313

Acetaldehyde is subject to reporting requirements of Section 313, Title III of SARA. 0.1 % de minimis concentration

Ethylene Oxide is subject to reporting requirements of Section 313, Title III of SARA. 0.1 % de minimis concentration

1,4-Dioxane is subject to reporting requirements of Section 313, Title III of SARA. 0.1 % de minimis concentration

Sodium Azide is subject to reporting requirements of Section 313, Title III of SARA. 1.0 % de minimis concentration

CERCLA RG's, 40 CFR 302.4

Acetaldehyde is listed.

Ethylene Oxide is listed.

1,4-Dioxane is listed.

Sodium Azide is listed.

California Proposition 65

Acetaldehyde has been identified by the State of California to cause cancer. The State of California has adopted a regulation which requires a warning be given to individual who may be exposed to chemicals identified by the State to cause cancer or reproductive harm.

WARNING: This product contains a chemical known to the State of California to cause cancer.

Ethylene Oxide has been identified by the State of California to cause cancer and reproductive harm. The State of California has adopted a regulation which requires a warning be given to individual who may be exposed to chemicals identified by the State to cause cancer or reproductive harm.

WARNING: This product contains a chemical known to the State of California to cause cancer and reproductive harm.

1,4-Dioxane has been identified by the State of California to cause cancer. The State of California has adopted a regulation which requires a warning be given to individual who may be exposed to chemicals identified by the State to cause cancer or reproductive harm.

WARNING: This product contains a chemical known to the State of California to cause cancer.

Section 15 Regulatory Information (Continued)

Massachusetts MSL

Acetaldehyde is listed.
Ethylene Oxide is listed.
1,4-Dioxane is listed.
Sodium Azide is listed.

New Jersey Dept. of Health RTK List

Acetaldehyde is listed.
Ethylene Oxide is listed.
1,4-Dioxane is listed.
Sodium Azide is listed.

Pennsylvania RTK

Acetaldehyde is listed.
Ethylene Oxide is listed.
1,4-Dioxane is listed.
Sodium Azide is listed.

EU Regulations

This SDS complies with EC Regulations 1907/2006 (REACH) and amendments.

Water Hazard Class (Germany) WGK 1, low water endangering

REACH 1907/2006 EC - Annex XIV - list of substances subject to authorization.

No ingredients listed.

According to EC Directives (1999/45/EC and 67/548 EEC)

Irritant

Xi



Risk and Safety Phrases

R43 May cause sensitization by skin contact.

S37 Wear suitable gloves.

S28 After contact with skin, wash immediately with plenty of soap and water.

Canada

This product is exempt from WHMIS label and SDS requirements.

PIN

Not applicable

Ingredients on Ingredient Disclosure List

Acetaldehyde
Ethylene Oxide
1,4-Dioxane
Sodium Azide

Ingredients with unknown toxicological properties

Product is exempt

15.2 Chemical Safety Assessment A Chemical Safety Assessment has not been carried out.

Some hazardous ingredients listed in Section 15 are below OSHAs and WHMIS' 1.0% w/w (0.1% for carcinogens) or EU's ingredient specific concentrations required for reporting in Section 3.

Section 16 Other Information

Beckman Coulter Safety Rating	Flammability: 0 Health: 2 Reactivity with Water: 0 Contact: 2	Code 0=None 1=Slight 2=Caution 3=Severe
--------------------------------------	--	---

Revision Changes Updated to GHS.

Hazard Class, hazard statements and risk phrase description from section 3

N - Dangerous for the environment
 T - Toxic
 T+ - Very toxic
 R23/24/25 Toxic by inhalation, in contact with skin and if swallowed.
 R34 Causes burns.
 R43 May cause sensitization by skin contact.
 R28 Very toxic if swallowed.
 R32 Contact with acids liberates very toxic gas.
 R50/53 Very toxic to aquatic organisms, may cause long-term adverse effects in the aquatic environment.
 Aquatic Acute 1 - Aquatic Hazard Acute, Category 1
 Acute Tox. Dermal 3 - Acute Toxicity Dermal, Category 3
 Acute Tox. Inhal. 3 - Acute Toxicity Inhalation, Category 3
 Acute Tox. Oral 2 - Acute Toxicity Oral, Category 2
 Acute Tox. Oral 3 - Acute Toxicity Oral, Category 3
 Aquatic Longterm 1 - Aquatic Hazard Long term, Category 1
 Skin Corr. 1B - Skin Corrosion Category 1B
 Skin Sens. 1 - Skin Sensitization Category 1
 H300 - Fatal if swallowed.
 H301 - Toxic if swallowed.
 H311 - Toxic in contact with skin.
 H314 - Causes severe skin burns and eye damage.
 H317 - May cause an allergic skin reaction.
 H331 - Toxic if inhaled.
 H400 - Very toxic to aquatic life.
 H410 - Very toxic to aquatic life with long lasting effects.

Abbreviations and Acronyms

ACGIH - American Conference of Governmental Industrial Hygienists
 ADR - European Agreement Concerning The International Carriage Of Dangerous Goods By Road
 CERCLA - The Comprehensive Environmental Response, Compensation, and Liability Act
 CLP - Classification, Labeling and Packaging
 DFGMAK - Republic Germany's maximum exposure limit
 GHS - Globally Harmonized System
 HCS - Hazard Communication Standard
 IARC - International Agency for Research on Cancer
 IATA - International Air Transport Association

Section 16 Other Information (Continued)

ICAO - International Civil Aviation Organization
IMDG - International Maritime Dangerous Goods
IOELVs - European Unions' Indicative Occupational Exposure Limit Values
NIOSH - National Institute for Occupational Safety and Health
NTP - National Toxicology Program
OSHA - Occupational Safety and Health Administration
PBT - Persistent bioaccumulative and toxic substances
SARA - Superfund Amendments and Reauthorization Act
TDG - Canadian Transportation Of Dangerous Goods Regulations.
UN GHS - United Nations Globally Harmonized System
US DOT - United States Department of Transportation
WHMIS - Workplace Hazardous Material Information System
vPvB - Very persistent and very bioaccumulative substances
LD50 - Lethal Dose, 50%
LC50 - Lethal Concentration, 50%

Access is a trademark of Beckman Coulter, Inc.

For further information, please contact your local Beckman Coulter, Inc. representative.

WHILE BECKMAN COULTER, INC. BELIEVES THE INFORMATION CONTAINED HEREIN IS VALID AND ACCURATE, MAKES NO WARRANTY OR REPRESENTATION AS TO ITS VALIDITY, ACCURACY, OR CURRENCY. BECKMAN COULTER, INC. SHALL NOT BE LIABLE OR OTHERWISE RESPONSIBLE IN ANY WAY FOR USE OF EITHER THIS INFORMATION OR MATERIALS TO WHICH IT APPLIES. DISPOSAL OF HAZARDOUS MATERIALS MAY BE SUBJECT TO LOCAL LAWS OR REGULATIONS.