

Safety Data Sheet (SDS)

Section 1: Chemical Product and Company Identification

Cat# : 4101, 4102, 4103, 4104

Part Name: Conflikt

Supplier: Decon Laboratories Inc.
460 Glennie Circle King of Prussia, Pa 19406
SDS Telephone # (610) 755-0800

Emergency Telephone Numbers

US Chemtrec: (800) 424-9300
Canada: (703) 527-3887

Identified uses: Laboratory use

Email Contact: cveloski@deconlabs.com

Section 2: Hazards Identification:

Hazard Overview

Physical hazards	Not classified	
Health hazards	Serious eye damage/eye irritation	Category 2B
Environmental hazards	Hazardous to the aquatic environment, acute hazard.	Category 2
	Hazardous to the aquatic environment, long-term hazard	Category 3
OSHA defined hazards	Not classified	

Label elements

Hazard symbol	None
Signal word	Warning
Hazard statement	Causes eye irritation. Toxic to aquatic life. Harmful to aquatic life with long lasting effects
Prevention	Wash thoroughly after handling. Avoid release to the environment.
Response	If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention.
Storage	Store away from incompatible materials.
Disposal	Dispose of contents/container in accordance with local/regional/national/international regulations.
Hazard(s) not otherwise classified (HNOC)	None known
Supplemental information	Not applicable.

Safety Data Sheet (SDS)

NFPA Rating

Hazard Ratings:

These ratings are Decon Laboratories Inc.'s own assessments of the properties of the material using the ANSI/NFPA 704 Standard. Additional information can be found by consulting in the NFPA published ratings lists (List 325 and list 49).

If no data is listed the information is not available

Health 1 Flammability 1 Reactivity 0

Section 3: Composition/ Information on ingredients

Mixtures

Chemical name	Common name and synonyms	CAS number	%
Diethylene glycol monobutyl ether		112-34-5	8
Tetra sodium ethylenediamine tetra acetic acid (Na4 EDTA)		64-02-8	1.6
Alkyl (68% C12, 32% C14) dimethyl ethylbenzyl ammonium chloride (Alternate CAS 68956-79-6)		85409-23-0	0.105
_Alkyl dimethyl benzyl ammonium chloride (C12-18)		68391-01-5	0.105
Other components below reportable levels			90.19

Section 4: First Aid Measures

Inhalation	Move to fresh air. Call a physician if symptoms develop or persist.
Skin contact	Wash off with soap and water. Get medical attention if irritation develops and persists.
Eye contact	Immediately flush eyes with plenty of water for at least 15 minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention.
Ingestion	Rinse mouth. Get medical attention if symptoms occur.
Most important symptoms/ effects, acute and delayed	Exposed individuals may experience eye tearing, redness, and discomfort.
Indication of immediate medical attention and special treatment needed.	If the product is ingested, probable mucosal damage may contraindicate the use of gastric lavage. Treat the affected person appropriately.
General information	Ensure that medical personnel are aware of the material(s) involved, and take precautions to protect themselves.

Section 5: Fire-Fighting Measures

Suitable extinguishing media Water fog.Foam. **Dry chemical powder. Carbon Dioxide (CO₂)**

Safety Data Sheet (SDS)

Unsuitable extinguishing media Do not use water jet as an extinguisher, as this will spread the fire.

Specific hazards arising from the chemical During fire, gases hazardous to health may be formed.

Special protective equipment and precautions for firefighters Self-contained breathing apparatus and full protective clothing must be worn in case of fire.

Fire-fighting equipment/instructions Move containers from fire area if you can do so without risk.

Specific methods Use standard firefighting procedures and consider the hazards of other involved materials.

General fire hazards No unusual fire or explosion hazards noted.

Section 6: Accidental Release measures

Personal precautions, protective equipment and emergency procedures Keep unnecessary personnel away. Keep people away from and upwind of spill/leak. Wear appropriate protective equipment and clothing during clean-up. Do not touch damaged containers or spilled material unless wearing appropriate protective clothing. Ensure adequate ventilation. Local authorities should be advised if significant spillages cannot be contained. For personal protection, see section 8 of the SDS.

Methods and materials for containment and cleaning up Large Spills: Stop the flow of material, if this is without risk. Dike the spilled material, where this is possible. Cover with plastic sheet to prevent spreading. Absorb in vermiculite, dry sand or earth and place into containers. Prevent product from entering drains. Following product recovery, flush area with water.

Small Spills: Wipe up with absorbent material (e.g. cloth, fleece). Clean surface thoroughly to remove residual contamination.

Never return spills to original containers for re-use. For waste disposal, see section 13 of the SDS.

Do not discharge effluent containing this product into lakes, streams, ponds, estuaries, oceans or other waters unless in accordance with the requirements of a National Pollutant Discharge Elimination System (NPDES) permit and the permitting authority has been notified in writing prior to discharge. Do not discharge effluent containing this product to sewer systems without previously notifying the local sewage treatment plant authority. For guidance, contact your State Water Board or Regional Office of the EPA.

Environmental precautions Avoid release to the environment. Contact local authorities in case of spillage to drain/aquatic environment. Prevent further leakage or spillage if safe to do so. Do not contaminate water. Avoid discharge into drains, water courses or onto the ground.

Safety Data Sheet (SDS)

Section 7: Handling and Storage

Precautions for safe handling Avoid contact with eyes. Avoid prolonged exposure. Provide adequate ventilation. Wear appropriate personal protective equipment. Observe good industrial hygiene practices. Avoid release to the environment. Do not empty into drains.

Conditions for safe storage, including any incompatibilities DO NOT CONTAMINATE WATER, FOOD OR FEED BY STORAGE OR DISPOSAL. PESTICIDE

STORAGE: Store in a dry place no lower in temperature than 50°F or higher than 120°F. Store in original tightly closed container. Store away from incompatible materials (see Section 10 of the SDS)

Section 8: Exposure Controls/ Personal Protection

Occupational exposure limits

US. ACGIH Threshold Limit Values

Components	Type	Value	Form
Diethylene glycol monobutyl ether (CAS 112-34-5)	TWA	10 ppm	Inhalable fraction and vapor.

Biological limit values No biological exposure limits noted for the ingredient(s).

Appropriate engineering controls Ensure adequate ventilation, especially in confined areas. Eye wash facilities and emergency shower must be available when handling this product.

Individual protection measures, such as personal protective equipment

Eye/face protection Wear safety glasses with side shields (or goggles).

Hand protection Wear protective gloves.

Skin protection

Other Wear suitable protective clothing

Respiratory protection In case of insufficient ventilation, wear suitable respiratory equipment

Thermal hazards Wear appropriate thermal protective clothing, when necessary.

General hygiene considerations Always observe good personal hygiene measures, such as washing after handling the material and before eating, drinking, and/or smoking. Routinely wash work clothing and protective equipment to remove contaminants.

Section 9: Physical and Chemical Properties

Appearance

Physical state	Liquid
Form	Liquid.
Color	Not available
Odor	Not available

Safety Data Sheet (SDS)

Odor threshold	Not available
pH	11.7000
Melting point/freezing point	> 32 °F (> 0 °C)
Initial boiling point and boiling range	> 212 °F (> 100 °C)
Flash point	> 201 °F (> 93.9 °C) Pensky-Martens Closed Cup
Evaporation rate	Not available
Flammability (solid, gas)	Not available.

Upper/lower flammability or explosive limits

Flammability limit – lower (%)	Not available.
Flammability limit - upper (%)	Not available
Explosive limit - lower (%)	Not available
Explosive limit - upper (%)	Not available

Vapor pressure	Not available.
Vapor density	Not available
Relative density	Not available

Solubility(ies)

Solubility (water)	Complete
---------------------------	----------

Auto-ignition temperature	Not available
Decomposition temperature	Not available
Viscosity	Not available

Other information

Density	8.38 lbs/gal @ 25 C
----------------	---------------------

Section 10: Stability and Reactivity:

Reactivity The product is stable and non-reactive under normal conditions of use, storage and transport

Chemical stability Material is stable under normal conditions

Possibility of hazardous reactions No dangerous reaction known under conditions of normal use.

Conditions to avoid Avoid temperatures exceeding the flash point. Contact with incompatible materials.

Incompatible materials Strong acids, alkalis and oxidizing agents.

Hazardous decomposition products Upon decomposition, this product may yield oxides of nitrogen and ammonia, carbon dioxide, carbon monoxide and other low molecular weight hydrocarbons.

Section 11: Toxicological Information

Safety Data Sheet (SDS)

Information on likely routes of exposure

Ingestion	Expected to be a low ingestion hazard.
Inhalation	Prolonged inhalation may be harmful
Skin contact	Causes mild skin irritation
Eye contact	Causes eye irritation.

Symptoms related to the physical, chemical and toxicological characteristics

Exposed individuals may experience eye tearing, redness, and discomfort.

Information on toxicological effects

Acute toxicity

Product

Detergent Disinfectant Pump Spray

Species

Test Results

Acute

Dermal

LD50

Rabbit

> 5 g/kg

Inhalation

LC50

Rat

> 2.43 mg/l

Oral

LD50

Rat

> 5 g/kg

Skin corrosion/irritation

Causes mild skin irritation

Serious eye damage/eye irritation

Causes eye irritation

Respiratory or skin sensitization

Respiratory sensitization

Not available

Skin sensitization

This product is not expected to cause skin sensitization

Germ cell mutagenicity

No data available to indicate product or any components present at greater than 0.1% are mutagenic or genotoxic.

Carcinogenicity

This product is not considered to be a carcinogen by IARC, ACGIH, NTP, or OSHA

US. OSHA Specifically Regulated Substances (29 CFR 1910.1001-1050)

Reproductive toxicity

This product is not expected to cause reproductive or developmental effects

Specific target organ toxicity - single exposure

Not applicable

Specific target organ toxicity - repeated exposure

Not applicable

Aspiration hazard

Not applicable

Safety Data Sheet (SDS)

Section 12: Ecological Information

Ecotoxicity Toxic to aquatic life. Harmful to aquatic life with long lasting effects

Components	Species	Test Results
Tetra sodium ethylenediamine tetra acetic acid (Na4 EDTA) (CAS 64-02-8)		
Aquatic		
Fish	LC50	Bluegill (Lepomis macrochirus) 472 - 500 mg/l, 96 hours

* Estimates for product may be based on additional component data not shown.

Persistence and degradability Expected to be readily biodegradable.

Bioaccumulative potential No data available.

Partition coefficient n-octanol / water (log Kow)

_Diethylene glycol monobutyl ether 0.56

Mobility in soil No data available.

Other adverse effects No other adverse environmental effects (e.g. ozone depletion, photochemical ozone creation potential, endocrine disruption, global warming potential) are expected from this component.

Section 13: Disposal Considerations

Disposal instructions PESTICIDE DISPOSAL - Pesticide wastes are acutely hazardous. Improper disposal of excess pesticide, spray mixture, or rinsate is a violation of Federal Law. If these wastes cannot be disposed of by use according to label instructions contact your State Pesticide or Environmental Control Agency, or the Hazardous Waste representative at the nearest EPA Regional Office for guidance.

CONTAINER DISPOSAL - Nonrefillable container. Do not reuse or refill container. Clean container promptly after emptying.

(For containers 5 gallons or less):

Triple rinse as follows: Fill container 1/4 full with water and recap. Agitate vigorously. Drain for 10 seconds after the flow begins to drip. Follow Pesticide Disposal instructions for rinsate disposal. Repeat procedure two more times. Then offer for recycling or reconditioning. If not available, puncture and dispose in a sanitary landfill.

(For containers greater than 5 gallons):

Triple rinse as follows: Fill container 1/4 full with water. Tip container on its side and roll it back and forth, ensuring at least one complete revolution for 30 seconds. Stand the container on its end and tip it back and forth several times. Turn the container over onto its other end and tip it back and forth several times. Follow Pesticide Disposal instructions for rinsate disposal. Repeat procedure two more times. Then offer for recycling or reconditioning. If not available, puncture and dispose in a sanitary landfill.

Hazardous waste code The waste code should be assigned in discussion between the user, the producer and the waste disposal company.

Safety Data Sheet (SDS)

Waste from residues / unused products Dispose of in accordance with local regulations. Empty containers or liners may retain some product residues. This material and its container must be disposed of in a safe manner (see: Disposal instructions).

Contaminated packaging Empty containers should be taken to an approved waste handling site for recycling or disposal. Since emptied containers may retain product residue, follow label warnings even after container is emptied.

Section 14: Transportation Information

General Note - Shipper/offeror is responsible for confirming appropriate proper shipping name, hazard classification(s), packing group, marking, labeling, placarding, packaging, and other regulatory requirements applicable to packages offered for transport. These requirements may vary depending on the quantity of material, packaging, mode of transportation, and carrier.

DOT Not regulated as dangerous goods.

IATA Not regulated as dangerous goods.

IMDG Not regulated as dangerous goods

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code Not available

Section 15: Regulatory Information

US federal regulations This product is a U.S. EPA registered pesticide.

This product is a "Hazardous Chemical" as defined by the OSHA Hazard Communication Standard, 29 CFR 1910.1200.

CERCLA Hazardous Substance List (40 CFR 302.4)

Not listed.

US. OSHA Specifically Regulated Substances (29 CFR 1910.1001-1050)

Not listed

Superfund Amendments and Reauthorization Act of 1986 (SARA

Hazard categories

Immediate Hazard – Yes

Delayed Hazard – No

Fire Hazard – No

Pressure Hazard – No

Reactivity Hazard – No

SARA 302 Extremely hazardous substance No

SARA 311/312 Hazardous chemical Yes
SARA 313 (TRI reporting)

Chemical name	CAS number	% by wt.
Acetaldehyde	75-07-0	<0.1

Other federal regulations

Clean Air Act (CAA) Section 112 Hazardous Air Pollutants (HAPs) List

Safety Data Sheet (SDS)

Not regulated
Clean Air Act (CCA) Section 112 (r) Accidental Release Prevention (40 CFR 68.130)
Not regulate
Safe Drinking Water Act (SDWA) Not regulated

FIFRA Information This chemical is a pesticide product registered by the Environmental Protection Agency and is subject to certain labeling requirements under federal pesticide law. These requirements differ from the classification criteria and hazard information required for safety data sheets, and for workplace labels of non-pesticide chemicals. Following is the hazard information as required on the pesticide label.

Signal word CAUTION
Hazard statement Causes moderate eye irritation

This pesticide is toxic to fish.

US state regulations

US. Massachusetts RTK - Substance List

Not regulated.

US. Rhode Island RTK

Not regulated

US. California Proposition 65

Subregistrant may elect not to use this warning if it conducts an assessment which determines such ingredient falls below the threshold action level requiring warning under California Proposition 65. Subregistrant does so at its own risk and liability, and hereby agrees to indemnify and hold Stepan harmless for any claim, cause of action, or damages arising from such decision not to warn.

WARNING: This product contains a chemical known to the State of California to cause cancer.

US - California Proposition 65 - CRT: Listed date/Carcinogenic substance

Acetaldehyde (CAS 75-07-0)

Listed: April 1, 1988

International Inventories

Country(s) or region	Inventory name	On inventory (yes/no)*
Australia	Australian Inventory of Chemical Substances (AICS)	No
Canada	Domestic Substances List (DSL)	No
Canada	Non-Domestic Substances List (NDSL)	No
China	Inventory of Existing Chemical Substances in China (IECSC)	No
Europe	European Inventory of Existing Commercial Chemical Substances (EINECS)	No
Europe	European List of Notified Chemical Substances (ELINCS)	No
Japan	Inventory of Existing and New Chemical Substances (ENCs)	No
Korea	Existing Chemicals List (ECL)	No
New Zealand	New Zealand Inventory	No

Safety Data Sheet (SDS)

Philippines	Philippine Inventory of Chemicals and Chemical Substances (PICCS)	No
United States & Puerto Rico	Toxic Substances Control Act (TSCA) Inventory	No

*A "Yes" indicates that all components of this product comply with the inventory requirements administered by the governing country(s)

A "No" indicates that one or more components of the product are not listed or exempt from listing on the inventory administered by the governing country(s).

Section 16: Other Information

Date of Issue: 1/1/1997

Date of Revision: 05/08/2015

Decon Laboratories, Inc. provides the information contained herein in good faith but makes no representation as to its comprehensiveness or accuracy. Individuals receiving this information must exercise their independent judgment in determining its appropriateness for a particular purpose. Decon Laboratories, Inc. makes no representations or warranties, either expressed or implied of merchantability, fitness for particular purposes with respect to the information set forth herein or to which the information refers. Accordingly, Decon Laboratories, Inc. will not be responsible for damages resulting from the use of or reliance upon this information.

End of Safety Data Sheet