

Safety Data Sheet

according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 12.16.2014

Page 1 of 7

Iodine, Gram Stain Solution

SECTION 1 : Identification of the substance/mixture and of the supplier

Product name : Iodine, Gram Stain Solution

Manufacturer/Supplier Trade name:

Manufacturer/Supplier Article number: S25941

Recommended uses of the product and restrictions on use:

Manufacturer Details:

AquaPhoenix Scientific, Inc
9 Barnhart Drive, Hanover, PA 17331
(717) 632-1291

Supplier Details:

Fisher Science Education
6771 Silver Crest Road, Nazareth, PA 18064
(724)517-1954

Emergency telephone number:

Fisher Science Education Emergency Telephone No.: 800-535-5053

SECTION 2 : Hazards identification

Classification of the substance or mixture:



Irritant

Skin irritation, category 2
Eye irritation, category 2A

Skin Irritation, Category 2

Eye Irritation, Category 2

Signal word :Warning

Hazard statements:

Causes serious eye irritation

Causes skin irritation

Precautionary statements:

If medical advice is needed, have product container or label at hand

Keep out of reach of children

Read label before use

Wash skin thoroughly after handling

Wear protective gloves/protective clothing/eye protection/face protection

IF ON SKIN: Wash with soap and water

Specific treatment (see supplemental first aid instructions on this label)

If skin irritation occurs: Get medical advice/attention

Take off contaminated clothing and wash before reuse

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do.

Continue rinsing

If eye irritation persists get medical advice/attention

Other Non-GHS Classification:

WHMIS

Safety Data Sheet


according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 12.16.2014

Page 2 of 7

Iodine, Gram Stain Solution

NFPA/HMIS



NFPA SCALE (0-4)

| | |
|---------------------|---|
| Health | 1 |
| Flammability | 0 |
| Physical Hazard | 0 |
| Personal Protection | X |

HMIS RATINGS (0-4)

SECTION 3 : Composition/information on ingredients

| Ingredients: | | |
|---------------------------|------------------|--------|
| CAS 7681-11-0 | Potassium Iodide | 3.05 % |
| CAS 7732-18-5 | DI Water | 95.1 % |
| CAS 7553-56-2 | Iodine | 1.85 % |
| Percentages are by weight | | |

SECTION 4 : First aid measures

Description of first aid measures

After inhalation: Move exposed individual to fresh air. Loosen clothing as necessary and position individual in a comfortable position. Seek medical advice if discomfort or irritation persists. If breathing difficult, give oxygen.

After skin contact: Wash affected area with soap and water. Rinse thoroughly. Seek medical attention if irritation persists or if concerned.

After eye contact: Protect unexposed eye. Rinse/flush exposed eye(s) gently using water for 15-20 minutes. Remove contact lens(es) if able to do so during rinsing. Seek medical attention if irritation persists or if concerned.

After swallowing: Rinse mouth thoroughly. Do not induce vomiting. Have exposed individual drink sips of water. Seek medical attention if irritation, discomfort or vomiting persists.

Most important symptoms and effects, both acute and delayed:

Irritation, Nausea, Headache, Shortness of breath.;

Indication of any immediate medical attention and special treatment needed:

If seeking medical attention, provide SDS document to physician.

SECTION 5 : Firefighting measures

Extinguishing media

Suitable extinguishing agents: If in laboratory setting, follow laboratory fire suppression procedures. Use appropriate fire suppression agents for adjacent combustible materials or sources of ignition. Water fog. Foam. Dry chemical powder. Carbon dioxide (CO₂).

For safety reasons unsuitable extinguishing agents:

Special hazards arising from the substance or mixture:

Combustion products may include carbon oxides or other toxic vapors. Thermal decomposition can lead to release of irritating gases and vapors.

Safety Data Sheet

according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 12.16.2014

Page 3 of 7

Iodine, Gram Stain Solution

Advice for firefighters:

Protective equipment: Use NIOSH-approved respiratory protection/breathing apparatus.

Additional information (precautions): Move product containers away from fire or keep cool with water spray as a protective measure, where feasible.

SECTION 6 : Accidental release measures

Personal precautions, protective equipment and emergency procedures:

Wear protective equipment. Use respiratory protective device against the effects of fumes/dust/aerosol. Keep unprotected persons away. Ensure adequate ventilation.

Environmental precautions:

Prevent from reaching drains, sewer or waterway. Collect contaminated soil for characterization per Section 13. Small quantities may be flushed to drains with plenty of water.

Methods and material for containment and cleaning up:

If in a laboratory setting, follow Chemical Hygiene Plan procedures. Collect liquids using vacuum or by use of absorbents. Place into properly labeled containers for recovery or disposal. If necessary, use trained response staff/contractor.

Reference to other sections:

SECTION 7 : Handling and storage

Precautions for safe handling:

Wash hands after handling. Follow good hygiene procedures when handling chemical materials. Use only in well ventilated areas. Avoid contact with eyes, skin, and clothing.

Conditions for safe storage, including any incompatibilities:

Provide ventilation for containers. Avoid storage near extreme heat, ignition sources or open flame. Store away from foodstuffs. Store away from oxidizing agents. Store in cool, dry conditions in well sealed containers. Keep container tightly sealed. Protect from freezing and physical damage.

SECTION 8 : Exposure controls/personal protection



Control Parameters:

7681-11-0, Potassium Iodide, ACS, ACGIH NIOSH 0.01 mg/m³

Appropriate Engineering controls:

Emergency eye wash fountains and safety showers should be available in the immediate vicinity of use/handling. Provide exhaust ventilation or other engineering controls to keep the airborne concentrations of vapor or dusts (total/respirable) below the applicable workplace exposure limits (Occupational Exposure Limits-OELs) indicated above. Use under a fume hood

Respiratory protection:

Use suitable respiratory protective device when high concentrations are present. For spills, respiratory protection may be advisable. Normal ventilation is adequate.

Protection of skin:

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation being used/handled. Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation.

Eye protection:

Safety glasses with side shields or goggles.

Safety Data Sheet

according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 12.16.2014

Page 4 of 7

Iodine, Gram Stain Solution

General hygienic measures: The usual precautionary measures are to be adhered to when handling chemicals. Wash hands before breaks and at the end of work. Do not inhale gases/fumes/dust/mist/vapor/aerosols. Avoid contact with the eyes and skin.

SECTION 9 : Physical and chemical properties

| | | | |
|---|-------------------------|--|---|
| Appearance (physical state,color): | Clear, colorless liquid | Explosion limit lower: Explosion limit upper: | Not determined Not determined |
| Odor: | Odorless | Vapor pressure: | Not Determined |
| Odor threshold: | Not determined | Vapor density: | Not determined |
| pH-value: | Not Determined | Relative density: | Approx 1 |
| Melting/Freezing point: | Approx 0°C | Solubilities: | Soluble in water |
| Boiling point/Boiling range: | Approx 100°C | Partition coefficient (n-octanol/water): | Not determined |
| Flash point (closed cup): | Not Determined | Auto/Self-ignition temperature: | Not determined |
| Evaporation rate: | Not determined | Decomposition temperature: | Not determined |
| Flammability (solid,gaseous): | Not determined | Viscosity: | a. Kinematic:Not determined b. Dynamic: Not determined |
| Density: Not determined | | | |

SECTION 10 : Stability and reactivity

Reactivity:Nonreactive under normal conditions.

Chemical stability:No decomposition if used and stored according to specifications.

Possible hazardous reactions:None under normal processing

Conditions to avoid:exposure to light.Incompatible Materials.

Incompatible materials:Strong acids.Strong bases.Strong oxidizers

Hazardous decomposition products:Hydrogen iodide. Iodine gas. May include oxides of iodine

SECTION 11 : Toxicological information

| | | |
|---|------------------------------|----------------------------|
| Acute Toxicity: | | |
| Oral: | POTASSIUM IODIDE (7681-11-0) | LD50 Rat: 285 mg/kg |
| Chronic Toxicity: No additional information. | | |
| Corrosion Irritation: | | |
| Dermal: | 7681-11-0 | Rabbit: causes irritation |
| Ocular: | 7681-11-0 | Rabbit: causes irritation |
| Sensitization: | | No additional information. |
| Single Target Organ (STOT): | | No additional information. |

Safety Data Sheet

according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 12.16.2014

Page 5 of 7

Iodine, Gram Stain Solution

| | |
|-------------------------------|----------------------------|
| Numerical Measures: | No additional information. |
| Carcinogenicity: | No additional information. |
| Mutagenicity: | No additional information. |
| Reproductive Toxicity: | No additional information. |

SECTION 12 : Ecological information

Ecotoxicity

Crustacea LC50 Zebra mussel (*Dreissena polymorpha*) 220 - 313 mg/l, 24 hours: 7681-11-0

Fish LC50 - *Oncorhynchus mykiss* (rainbow trout) - 2,190 mg/l - 96 h: 7681-11-0

Persistence and degradability:

Bioaccumulative potential: Not Bioaccumulative.

Mobility in soil:

Other adverse effects:

SECTION 13 : Disposal considerations

Waste disposal recommendations:

Product/containers must not be disposed together with household garbage. It is the responsibility of the waste generator to properly characterize all waste materials according to applicable regulatory entities (US 40CFR262.11). Consult federal state/ provincial and local regulations regarding the proper disposal of waste material that may incorporate some amount of this product. Small amounts may be flushed with water to sewer. Larger volumes must be sent to approved plant for destruction.

SECTION 14 : Transport information

UN-Number

Not Regulated.

UN proper shipping name

Not Regulated.

Transport hazard class(es)

Packing group: Not Regulated

Environmental hazard:

Transport in bulk:

Special precautions for user:

SECTION 15 : Regulatory information

United States (USA)

SARA Section 311/312 (Specific toxic chemical listings):

Acute

SARA Section 313 (Specific toxic chemical listings):

None of the ingredients is listed

RCRA (hazardous waste code):

None of the ingredients is listed

TSCA (Toxic Substances Control Act):

Safety Data Sheet

according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 12.16.2014

Page 6 of 7

Iodine, Gram Stain Solution

All ingredients are listed.

CERCLA (Comprehensive Environmental Response, Compensation, and Liability Act):

None of the ingredients is listed

Proposition 65 (California):

Chemicals known to cause cancer:

None of the ingredients is listed

Chemicals known to cause reproductive toxicity for females:

None of the ingredients is listed

Chemicals known to cause reproductive toxicity for males:

None of the ingredients is listed

Chemicals known to cause developmental toxicity:

None of the ingredients is listed

Canada

Canadian Domestic Substances List (DSL):

All ingredients are listed.

Canadian NPRI Ingredient Disclosure list (limit 0.1%):

None of the ingredients is listed

Canadian NPRI Ingredient Disclosure list (limit 1%):

None of the ingredients is listed

SECTION 16 : Other information

This product has been classified in accordance with hazard criteria of the Controlled Products Regulations and the SDS contains all the information required by the Controlled Products Regulations. Note: The responsibility to provide a safe workplace remains with the user. The user should consider the health hazards and safety information contained herein as a guide and should take those precautions required in an individual operation to instruct employees and develop work practice procedures for a safe work environment. The information contained herein is, to the best of our knowledge and belief, accurate. However, since the conditions of handling and use are beyond our control, we make no guarantee of results, and assume no liability for damages incurred by the use of this material. It is the responsibility of the user to comply with all applicable laws and regulations applicable to this material.

GHS Full Text Phrases:

Abbreviations and acronyms:

IMDG: International Maritime Code for Dangerous Goods

PNEC: Predicted No-Effect Concentration (REACH)

CFR: Code of Federal Regulations (USA)

SARA: Superfund Amendments and Reauthorization Act (USA)

RCRA: Resource Conservation and Recovery Act (USA)

TSCA: Toxic Substances Control Act (USA)

NPRI: National Pollutant Release Inventory (Canada)

DOT: US Department of Transportation

IATA: International Air Transport Association

GHS: Globally Harmonized System of Classification and Labelling of Chemicals

ACGIH: American Conference of Governmental Industrial Hygienists

CAS: Chemical Abstracts Service (division of the American Chemical Society)

Safety Data Sheet

according to 29CFR1910/1200 and GHS Rev. 3

Effective date : 12.16.2014

Page 7 of 7

Iodine, Gram Stain Solution

NFPA: National Fire Protection Association (USA)

HMIS: Hazardous Materials Identification System (USA)

WHMIS: Workplace Hazardous Materials Information System (Canada)

DNEL: Derived No-Effect Level (REACH)

Effective date : 12.16.2014

Last updated : 05.30.2015