

GG&B COMPANY
3411 McNIEL AVENUE
SUITE 302
WICHITA FALLS, TX 76308

SAFETY DATA SHEET (SDS)

#1 - IDENTIFICATION

Product Identifier – Quick Slide Stain Pack

Product Code – 29000

Synonyms – Wright Giemsa Rinse

Recommended use - Professional, scientific and medical-technical activities

Recommended restrictions - None known

Manufacturer:

Company Name - ENG Scientific, Inc. (a.k.a. A.J.P Scientific, Inc.)

82 Industrial St East

Clifton, New Jersey USA 07012-1708

Telephone - Telephone (973) 472 - 7200

Toll Free (800) 922 - 0223

Fax (973) 472 - 9460

Emergency Contact Chemtrec (800) 424-9300

#2 - Hazard(s) Identification

Physical hazards	Flammable liquid	Category 1
	Health hazard	Category 0

Label elements



Signal word	Flammable, Poison
Hazard statement	Flammable liquid
Precautions	Keep away from heat/open flames/hot surfaces
Storage	Store at room temperature
Supplemental info	Avoid release to the environment
Disposal	Dispose of contents-according to local regulations

NFPA Rating	- HEALTH:	0
	- FLAMMABILITY:	1
	- REACTIVITY:	0

#3 – Composition

Chemical name: Methanol CAS # 67-56-1

#4 – First-aid measures

Inhalation: Remove to fresh air.

Eye contact: Remove contact lenses if worn. In case of contact, immediately flush eyes with plenty of clean running water for at least 15 minutes, lifting the upper and lower eyelids occasionally. Obtain medical attention.

Skin contact: Wash skin and clothing with soap and water.

Swallowing: Swallowing methanol is potentially life threatening. Onset of symptoms may be delayed for 18 to 24 hours. Do not induce vomiting.

#5 – Fire-fighting measures

Extinguishing media: Use water spray, carbon dioxide, dry chemical, alcohol-type, or universal type foams applied by manufacturer's recommended technique.

Special fire fighting procedures: Methanol burns with a clean clear flame that is almost invisible in daylight. Closed containers may rupture violently and suddenly release methanol when exposed to fire or excessive heat for a sufficient period of time. Vapors are slightly heavier than air and may travel long distances toward sources of ignition.

#6 – Accidental release measures

Steps to be taken in case material is released or spilled: Extinguish all sources of ignition; ventilate the area; flood or flush the area with water. Prevent spilled methanol from entering sewers, confined spaces, drains, or waterways. Put material in suitable, covered, labeled containers. Flush area with water.

#7 – Handling and storage

Handling: No smoking or open flame in storage, use or handling areas. Use explosion proof electrical equipment. Ensure proper grounding procedures are in place.

Storage: Store at room temperature and in a tightly closed container.

#8 – Exposure controls/personal protection

Respiratory protection: NOISH - approved respirator.

Ventilation: General (mechanical) room ventilation recommended. Special local ventilation also recommended.

Protective gloves: Impermeable.

Eye protection: Chemical splash goggles or face shield.

Other protective equipment: Lab coat. Eye bath, safety shower.

#9 – Physical and chemical properties

Appearance: Clear, colorless liquid.

Odor: Mild characteristic alcohol odor. **Boiling point:** 64.7°C

Specific gravity: (H₂O = 1) 0.791 @ 20°C

Vapor density: (air = 1) 1.105 **Solubility in H₂O** = Complete **Freezing point:** -97.8°C **Evaporation Rate** = 5.8

Flash point: (test methods) 108°F, Tag Closed Cup ASTM D 56

#10 – Stability and reactivity

Stability: Stable as supplied.

Incompatibility: Avoid contact with strong oxidizers, strong mineral or organic acids, and strong bases. May be corrosive to lead, aluminum, magnesium, and platinum.

Conditions to avoid: Avoid contact with sparks, heat, open flame, or ignition sources.

Hazardous polymerization: Will not occur.

#11 - Toxicology information

Route of Exposure - Inhalation

Prolonged inhalation of vapors causes dizziness, nausea, visual impairment, respiratory failure, muscular non-coordination and narcosis. Inhalation of high concentrations for prolonged periods has resulted in death. Liver damage has resulted from prolonged or repeated inhalation of vapors.

Route of Exposure - Skin

Essentially non-irritating. Repeated or prolonged contact causes drying, brittleness, cracking and irritation. Slightly toxic to animals by absorption.

First Aid - Skin

Remove contaminated clothing and wash contaminated skin with large amounts water. If irritation persists, contact a physician,

Route of Exposure - Eyes

May cause eye injury which may persist for several days. Liquid and vapor in high concentrations, causes irritation, tearing and burning sensation.

First Aid - Eyes

Flush eyes with water for at least 15 minutes. Contact a physician

Route of Exposure - Ingestion

Poisonous if swallowed. Can affect the optic nerve resulting in blindness. Can cause mental sluggishness, nausea and vomiting leading to severe illness, possibly death (in humans). Practically non-toxic to animals.

Health Conditions Aggravated By Exposure

Methanol exposure may aggravate existing eye, skin, kidney and liver disorders. Pre placement and annual medical examinations are recommended for workers who normally handle methanol with emphasis on neurological, visual, liver and kidney functions. Monitoring of air in the work place is recommended to maintain methanol vapors below recommended TLV.

#12 – Ecological information

DO NOT discharge into sewer or waterways. Small spills should be flushed with large quantities of water. Large spills should be collected for disposal. Flush with water, controlled burning.

#13 – Disposal considerations

Disposal of hazardous material may be subject to federal, state or local regulations. The information contained in this document is believed to be factual, however, nothing presented herein is to be taken as a warranty or representation for which ENG Scientific, Inc.(a.k.a. A.J.P. Scientific, Inc.) shall be responsible.

#14 – Transport information

PROPER SHIPPING NAME: N/A

UN NUMBER: N/A

HAZARD CLASS: 0

PACKING GROUP: N/A

DOT SHIPPING LABEL:

N/A For Limited Quantity

#15 – Regulatory information

US federal regulations All components are on the U.S. EPA TSCA Inventory List.

TSCA Section 12(b)Export Notification (40 CFR 707, SUBPT. D) Not regulated

US. OSHA Specifically Regulated Substances (29 CFR 1910.1001-1050) Not on regulatory list

CERCLA Hazardous Substance List (40 CFR 302.4) Methanol (CAS 67-56-1)

Superfund Amendments and Reauthorization Act of 1986 (SARA)

Hazard categories

Immediate Hazard – Yes

Delayed Hazard –Yes

Fire Hazard – Yes

Pressure Hazard – No

Reactivity Hazard – No

SARA 302 Extremely hazardous substance No

Other federal regulations

Clean Air Act (CAA) Section 112 Hazardous Air Pollutants (HAPs) List: Not regulated

Clean Air Act (CAA) Section 112® Accidental Release Prevention (40 CFR 68.130): Not regulated

Clean Water Act (CWA) Section 112® (40 CFR 68.130) : Hazardous substance

Safe Drinking Water Act (SDWA): Not regulated.

Drug enforcement Administration (DEA). List 2 Essential Chemicals (21 XFR 1310.02(b) and 1310.04(f)(2) and Chemical Code Number : Not listed.

Drug Enforcement Administration (DEA). List 1 & 2 Exempt Chemical Mixtures (21 CFR 1310.12(c) :
Not regulated.

DEA Exempt Chemical Mixtures Code Number: Not regulated.

Food and Drug Administration (FDA): N/A

US state regulation: California Safe Drinking Water and Toxic Enforcement Act of 1986 (Proposition 65):

This material is not known to contain and chemicals currently listed as carcinogens or reproductive toxins.

US. Massachusetts RTK – Substance List: Methanol (CAS 67-19-7)

US. New Jersey Worker and Community Right-to-Know Act: Not regulated

US. Pennsylvania RTK – Hazardous Substances: Methanol (CAS 67-19-7)

US. Rhode Island RTK: Methanol (CAS 67-19-7)

US California Proposition 65:

US-California Proposition 65 – Carcinogens and Reproductive Toxicity (CRT): Not listed.

International Inventories

Country/Region	Inventory name	On Inventory*
Canada	Non-domestic Substances List(NDSL)	Yes
Europe	European Inventory of Existing Chemical Substances (EINECS)	Yes
U.S.	Toxic Substances Control Act (TSCA) Inventory	Yes
Puerto Rico	Toxic Substances Control Act (TSCA) Inventory	Yes

** Yes complies with the inventory requirements administered by the governing country.*

#16 – Other information

Issue date March 2015

Version # 01

Further information The information contained in this document is believed to be factual, however, nothing presented herein is to be taken as a warranty or representation for which ENG Scientific, Inc. (a.k.a., A.J.P. Scientific, Inc.) shall be responsible.