

SAFETY DATA SHEET

Version 5.2
Revision Date 08/27/2014
Print Date 03/09/2015

1. PRODUCT AND COMPANY IDENTIFICATION**1.1 Product identifiers**

Product name : naphthol as-D chloroacetate (specific) esterase kit
for 90C2

Product Number : 90C2
Brand : Sigma

1.2 Relevant identified uses of the substance or mixture and uses advised against

Identified uses : Laboratory chemicals, Manufacture of substances

1.3 Details of the supplier of the safety data sheet

Company : Sigma-Aldrich
3050 Spruce Street
SAINT LOUIS MO 63103
USA

Telephone : +1 800-325-5832
Fax : +1 800-325-5052

1.4 Emergency telephone number

Emergency Phone # : (314) 776-6555

This is a summary MSDS for a kit, for the full MSDS for each of the components listed in section 16 please visit our website.

2. HAZARDS IDENTIFICATION**2.1 Classification of the substance or mixture****GHS Classification in accordance with 29 CFR 1910 (OSHA HCS)**

Acute toxicity, Oral (Category 4), H302
Acute toxicity, Inhalation (Category 4), H332
Acute toxicity, Dermal (Category 4), H312
Skin irritation (Category 2), H315
Serious eye damage (Category 1), H318
Respiratory sensitisation (Category 1), H334
Skin sensitisation (Category 1), H317
Carcinogenicity (Category 1B), H350
Reproductive toxicity (Category 1B), H360

For the full text of the H-Statements mentioned in this Section, see Section 16.

2.2 GHS Label elements, including precautionary statements

Pictogram



Signal word

Danger

Hazard statement(s)
H302 + H312 + H332
H315
H317

Harmful if swallowed, in contact with skin or if inhaled
Causes skin irritation.
May cause an allergic skin reaction.

H318	Causes serious eye damage.
H334	May cause allergy or asthma symptoms or breathing difficulties if inhaled.
H350	May cause cancer.
H360	May damage fertility or the unborn child.
Precautionary statement(s)	
P201	Obtain special instructions before use.
P202	Do not handle until all safety precautions have been read and understood.
P261	Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.
P264	Wash skin thoroughly after handling.
P270	Do not eat, drink or smoke when using this product.
P271	Use only outdoors or in a well-ventilated area.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves/ eye protection/ face protection.
P301 + P312	IF SWALLOWED: Call a POISON CENTER or doctor/ physician if you feel unwell.
P302 + P352	IF ON SKIN: Wash with plenty of soap and water.
P304 + P340	IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.
P305 + P351 + P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P310	Immediately call a POISON CENTER or doctor/ physician.
P322	Specific measures (see supplemental first aid instructions on this label).
P330	Rinse mouth.
P333 + P313	If skin irritation or rash occurs: Get medical advice/ attention.
P362	Take off contaminated clothing and wash before reuse.
P405	Store locked up.
P501	Dispose of contents/ container to an approved waste disposal plant.

2.3 Hazards not otherwise classified (HNOC) or not covered by GHS - none

3. COMPOSITION/INFORMATION ON INGREDIENTS

Refer to component MSDS

4. FIRST AID MEASURES

Refer to component MSDS

5. FIREFIGHTING MEASURES

Refer to component MSDS

6. ACCIDENTAL RELEASE MEASURES

Refer to component MSDS

7. HANDLING AND STORAGE

7.1 Precautions for safe handling

For precautions see section 2.2.

7.2 Conditions for safe storage, including any incompatibilities

Keep containers tightly closed in a dry, cool and well-ventilated place. Store in cool place.

Recommended storage temperature 2 - 8 °C

7.3 Specific end use(s)

Apart from the uses mentioned in section 1.2 no other specific uses are stipulated

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Refer to component MSDS

9. PHYSICAL AND CHEMICAL PROPERTIES

Refer to component MSDS

10. STABILITY AND REACTIVITY

Refer to component MSDS

11. TOXICOLOGICAL INFORMATION

Refer to component MSDS

12. ECOLOGICAL INFORMATION

Refer to component MSDS

13. DISPOSAL CONSIDERATIONS

Refer to component MSDS

14. TRANSPORT INFORMATION**DOT (US)**

UN number: 3316 Class: 9 Packing group: III
Proper shipping name: Chemical kits
Marine pollutant: No
Poison Inhalation Hazard: No

IMDG

UN number: 3316 Class: 9 Packing group: III EMS-No: F-A, S-P
Proper shipping name: CHEMICAL KIT
Marine pollutant: No

IATA

UN number: 3316 Class: 9 Packing group: III
Proper shipping name: Chemical kit

15. REGULATORY INFORMATION**SARA 311/312 Hazards**

Acute Health Hazard, Chronic Health Hazard

16. OTHER INFORMATION**Kit Components:**

Trizmal™ buffer	SIGMA	903C	Carc. 2; H351
N,N-Dimethylformamide	SIGMA	9010	Flam. Liq. 3; Acute Tox. 4; Eye Irrit. 2A; Repr. 1B; H226, H312 + H332, H319, H360
2-Methoxy-5-methyl-4-[(4-methyl-2-nitrophenyl)azo]benzenediazonium tetrachlorozincate (2:1)	SIGMA	9015	Acute Tox. 4; Carc. 1B; H302 + H312 + H332, H350
Acid hematoxylin solution	SIGMA	2852	Acute Tox. 4; Skin Irrit. 2; Eye Irrit. 2A; Resp. Sens. 1; Skin Sens. 1; H302, H315, H317, H319, H334
Citrate Concentrate	SIGMA	3861	Skin Irrit. 2; Eye Dam. 1; H315, H318

HMIS Rating

Health hazard:	3
Chronic Health Hazard:	*
Flammability:	0
Physical Hazard	1

NFPA Rating

Health hazard:	3
Fire Hazard:	0
Reactivity Hazard:	0

Further information

The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. The information in this document is based on the present state of our knowledge and is applicable to the product with regard to appropriate safety precautions. It does not represent any guarantee of the properties of the product. Sigma-Aldrich Corporation and its Affiliates shall not be held liable for any damage resulting from handling or from contact with the above product. See www.sigma-aldrich.com and/or the reverse side of invoice or packing slip for additional terms and conditions of sale.

Copyright 2014 Sigma-Aldrich Co. LLC. License granted to make unlimited paper copies for internal use only.

Preparation Information

Sigma-Aldrich Corporation
Product Safety – Americas Region
1-800-521-8956

Version: 5.2

Revision Date: 08/27/2014

Print Date: 03/09/2015