

SAFETY DATA SHEET

Version 5.4
Revision Date 08/21/2014
Print Date 03/09/2015

1. PRODUCT AND COMPANY IDENTIFICATION**1.1 Product identifiers**

Product name : Naphthol AS-D Chloroacetate (Specific Esterase) Kit
Product Number : 91C
Brand : Sigma

1.2 Relevant identified uses of the substance or mixture and uses advised against

Identified uses : Laboratory chemicals, Manufacture of substances

1.3 Details of the supplier of the safety data sheet

Company : Sigma-Aldrich
3050 Spruce Street
SAINT LOUIS MO 63103
USA
Telephone : +1 800-325-5832
Fax : +1 800-325-5052

1.4 Emergency telephone number

Emergency Phone # : (314) 776-6555

This is a summary MSDS for a kit, for the full MSDS for each of the components listed in section 16 please visit our website.

2. HAZARDS IDENTIFICATION**2.1 Classification of the substance or mixture****GHS Classification in accordance with 29 CFR 1910 (OSHA HCS)**

Skin corrosion (Category 1A), H314
Serious eye damage (Category 1), H318
Skin sensitisation (Category 1), H317

For the full text of the H-Statements mentioned in this Section, see Section 16.

2.2 GHS Label elements, including precautionary statements

Pictogram



Signal word

Danger

Hazard statement(s)

H314
H317

Causes severe skin burns and eye damage.
May cause an allergic skin reaction.

Precautionary statement(s)

P261
P264
P272
P280

Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.
Wash skin thoroughly after handling.
Contaminated work clothing should not be allowed out of the workplace.
Wear protective gloves/ protective clothing/ eye protection/ face protection.

P301 + P330 + P331

IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.

P303 + P361 + P353	IF ON SKIN (or hair): Remove/ Take off immediately all contaminated clothing. Rinse skin with water/ shower.
P304 + P340	IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.
P305 + P351 + P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P310	Immediately call a POISON CENTER or doctor/ physician.
P321	Specific treatment (see supplemental first aid instructions on this label).
P333 + P313	If skin irritation or rash occurs: Get medical advice/ attention.
P363	Wash contaminated clothing before reuse.
P405	Store locked up.
P501	Dispose of contents/ container to an approved waste disposal plant.

2.3 Hazards not otherwise classified (HNOC) or not covered by GHS - none

3. COMPOSITION/INFORMATION ON INGREDIENTS

Refer to component MSDS

4. FIRST AID MEASURES

Refer to component MSDS

5. FIREFIGHTING MEASURES

Refer to component MSDS

6. ACCIDENTAL RELEASE MEASURES

Refer to component MSDS

7. HANDLING AND STORAGE

7.1 Precautions for safe handling

For precautions see section 2.2.

7.2 Conditions for safe storage, including any incompatibilities

Keep container tightly closed in a dry and well-ventilated place.

Recommended storage temperature 2 - 8 °C

7.3 Specific end use(s)

Apart from the uses mentioned in section 1.2 no other specific uses are stipulated

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Refer to component MSDS

9. PHYSICAL AND CHEMICAL PROPERTIES

Refer to component MSDS

10. STABILITY AND REACTIVITY

Refer to component MSDS

11. TOXICOLOGICAL INFORMATION

Refer to component MSDS

12. ECOLOGICAL INFORMATION

Refer to component MSDS

13. DISPOSAL CONSIDERATIONS

Refer to component MSDS

14. TRANSPORT INFORMATION

DOT (US)

UN number: 3316 Class: 9 Packing group: II
Proper shipping name: Chemical kits
Marine pollutant: No
Poison Inhalation Hazard: No

IMDG

UN number: 3316 Class: 9 Packing group: II EMS-No: F-A, S-P
Proper shipping name: CHEMICAL KIT
Marine pollutant: No

IATA

UN number: 3316 Class: 9 Packing group: II
Proper shipping name: Chemical kit

15. REGULATORY INFORMATION

SARA 311/312 Hazards

Acute Health Hazard

16. OTHER INFORMATION

Kit Components:

Citrate Solution	SIGMA	915	
Fast Red Violet LB Base solution	SIGMA	912	Skin Corr. 1A; Eye Dam. 1; H314
Naphthol AS-D Chloroacetate solution	SIGMA	911	Skin Sens. 1; H317

HMIS Rating

Health hazard: 3
Chronic Health Hazard:
Flammability: 0
Physical Hazard 0

NFPA Rating

Health hazard: 3
Fire Hazard: 0
Reactivity Hazard: 0

Further information

Copyright 2014 Sigma-Aldrich Co. LLC. License granted to make unlimited paper copies for internal use only. The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. The information in this document is based on the present state of our knowledge and is applicable to the product with regard to appropriate safety precautions. It does not represent any guarantee of the properties of the product. Sigma-Aldrich Corporation and its Affiliates shall not be held liable for any damage resulting from handling or from contact with the above product. See www.sigma-aldrich.com and/or the reverse side of invoice or packing slip for additional terms and conditions of sale.

Preparation Information

Sigma-Aldrich Corporation
Product Safety – Americas Region
1-800-521-8956

Version: 5.4

Revision Date: 08/21/2014

Print Date: 03/09/2015