

Safety Data Sheet

Complies with EC no. 1907/2006
Issue Date: 11/15/2013

Section 1: Chemical Product and Company Identification

Page 1 of 4

Cat #: PG, PG-16

Part Name: PARAGARD PARAFFIN REPELLENT

Supplier: General Data Healthcare (www.general-data.com/healthcare)
4354 Ferguson Drive
Cincinnati, OH 45245 USA
Phone: 844.643.1129
24 hr. Emergency Contact 1.800.255.3924

Identified uses: Laboratory use, manufacture of substances

Section 2: Hazards Identification

Hazard Overview

Causes respiratory tract irritation.

Combustible liquid and vapor.

GHS Classification

Acute Toxicity Inhalation Category 5

Flammable Liquids Category 4

Signal word: Danger



Hazard and Precautionary Statements

- H227 Combustible liquid
H332 Harmful if inhaled.
P285 In case of inadequate ventilation wear respiratory protection.
P301A IF SWALLOWED do not induce vomiting. Do not give anything to drink. Obtain medical attention without delay.
P302+P352 IF ON SKIN: Wash with plenty of soap and water.
P304+P341 IF INHALED: If breathing is difficult, remove victim to fresh air and keep at rest in a position comfortable for breathing.
P314 Get medical advice/attention if you feel unwell.
P351 Rinse cautiously with water for several minutes.
P361 Remove/Take off immediately all contaminated clothing.

NFPA Rating

Hazard Ratings:

These ratings are private assessments of the properties of the material using the ANSI/NFPA 704 Standard. Additional information can be found by consulting in the NFPA published ratings lists (List 325 and List 49).

If no data is listed the information is not available.

Health	Flammability	Reactivity
1	2	0

Section 3: Composition/ Information on Ingredients

Note: Items listed with a CASRN... number have no CAS# available.

Item#	Name	EINECS	CAS#	% in product
1	Methyl salicylate	204-317-7	0000119368	0 - 5
2	Mineral oil	232-455-8	0008042475	51 - 60
3	Naphtha, petroleum, hydrotreated heavy	265-150-3	0064742489	41 - 50

Section 4: First Aid Measures

Chemical Stability: stable
Conditions to Avoid: flame
Incompatibility with other materials: none
Hazardous Decomposition Products: none
Hazardous Polymerization: will not occur

Section 11: Toxicological Information

Acute Data: no data
Subchronic data: no data

Section 12: Ecological Information

no data

Section 13: Disposal Considerations

The following chart lists the status of the chemical and its components in reference to 40 CFR Part 261.33. If the product is listed by code number the substance may be subject to special federal and state disposal regulations. If no codes are listed the material must be disposed in compliance with all Federal, State and Local Regulations.

CAS#	Waste Code	Regulated Name
0000119368	not listed	not listed
0008042475	not listed	not listed
0064742489	not listed	not listed

Section 14: Transportation Data

Proper Shipping Name: Combustible liquid, n.o.s
Chemical Name: PETROLEUM DISTILLATES

UN
Class
PG

Section 15: Regulatory Information

All components of this product are on the TSCA public inventory.

All components of this product are on the TSCA public inventory.

Prop 65 - Column A identifies those items which are known to the State of California to cause cancer. Column B identified items which are known to the State of California to cause reproductive toxicity.

CAS#	Column A	Column B
0000119368	no	no
0008042475	no	no
0064742489	no	no

State Regulatory Information :If a CAS# is listed below this material is subject to the listed state right-to-know requirements.

CAS#	
0000119368	not listed
0008042475	not listed
0064742489	not listed

SARA Toxic Release Chemicals(as defined in Section 313 of SARA Title III)

This list identifies the toxic chemicals, including their de minimis concentrations for which reporting is required under Section 313 of the Emergency Planning and Community Right-to-Know Act (EPCRA). The list is also referred to as the Toxics Release Inventory (TRI) List.

CAS#	Regulated name	de minimis conc. %	Rep. Thres.
0000119368	not listed	not listed	not listed
0008042475	not listed	not listed	not listed
0064742489	not listed	not listed	not listed

SARA Extremely Hazardous Substances and TPQs

This list includes hazardous chemicals as defined in 29 CFR 1910.1200(c); and extremely hazardous substances regulated under Section 302 of SARA Title III with their TPQs (in pounds), as listed in 40 CFR 355, Appendices A and B.

CAS#	Regulated name	TPQ (pounds)	EHS-RQ(pounds)
0000119368	not listed	not listed	not listed

0008042475	not listed	not listed	not listed
0064742489	not listed	not listed	not listed

CERCLA

The hazardous substances, and their reportable quantities (RQs) are listed in the federal regulations at 40 CFR Part 302, Table 302.4. Release of a CERCLA hazardous substance in an amount equal to or greater than its RQ, in any 24-hour period, must be reported to the National Response Center at (800) 424-8802.

CAS#	Regulated name	RQ (pounds)
0000119368	Not listed	Not listed
0008042475	Not listed	Not listed
0064742489	Not listed	Not listed

Section 16: Other Information

General Data Healthcare provides the information contained herein in good faith but makes no representation as to its comprehensiveness or accuracy. Individuals receiving this information must exercise their independent judgment in determining its appropriateness for a particular purpose. General Data Healthcare makes no representations or warranties, either expressed or implied of merchantability, fitness for particular purposes with respect to the information set forth herein or to which the information refers. Accordingly, General Data Healthcare will not be responsible for damages resulting from the use of or reliance upon this information.

END OF MSDS