

# Bio Lion<sup>®</sup> XC-L5

## CENTRIFUGE & ROTOR COMBO

This all purpose low speed centrifuge is quiet, compact, and convenient which makes it ideal for any lab. With swing bucket rotors that can accommodate up to 50mL tubes, this versatile centrifuge can be used in applications from clinical to biochemical for qualitative analysis of serum, plasma, blood, and urine.

### WEIGHT & DIMENSIONS

Max Speed: 5000r/min

Max RCF: 4390 X g

Max Capacity: 500mL X 4

Swing Rotor (CAS S0171B):  
50mL X 8 (4000r/min)

RCF: 2810 X g


Timer Range: 0~99min

Noise: <65dB

Power Supply: AC110V 60Hz 10A

Dimensions: 460 X 550 X 340mm (LxWxH)  
18" X 21.5" X 13.5"

Shipping Weight: 100 lbs



CAS XCL5  
CENTRIFUGE

CAS S0171B  
ROTOR

### FEATURES

- LED Screen
- Brushless Motor
- Low Maintenance
- Programmable
- Automatic RCF Display
- Comes with balance control, lid lock, and temperature/speed protection devices

**Mercedes**  
MEDICAL

PHONE: 800.331.2716 | FAX: 800.359.8807 | WEB: MERCEDESMEDICAL.COM

All offers are subject to availability and are subject to change without notice. Please contact one of our Customer Care Representatives at 800.331.2716 to confirm pricing. Mercedes Medical expressly disclaims any affiliation with or sponsorship referred to within. The names of various manufacturers, their instruments and their products referred to herein may be protected by trademark or other law, and are used herein solely for purpose of reference. All trademarks are property of their respective owners. To stop receiving faxes, please fill out the information below and fax back to the number above.

63487